

AlphaProTech

NIOSH Approved
N-95 Particulate Respirator
Part No. 695
Contents: 35
 *Fit Test required per NIOSH

N-95 Particulate Respirator
Positive Facial Lock®
 With Magic Arch®
 One Size Fits All*

Positive Facial Lock®
 With Magic Arch®

CUSTOM FIT MASK

Directions for Use:



1. Grasp nosepiece at top of mask with one hand. Slide finger of other hand behind chin piece.



2. Open mask by pulling down and out on chin piece until mask is fully open. It is IMPORTANT that mask is fully open.



3. Place chin in pocket created. Stretch bands over back of head, positioning bottom band behind neck below ears and top band near crown of head.



4. Conform nosepiece to nose and face contours making sure Magic Arch® is fully extended away from nose and mouth.



5. Holding mask at nosepiece, place one finger inside the mask on the malleable aluminum chin piece, pull down gently so that mask is snug against face.



6. Pinch any excess material together under chin and fold over or twist.



7. Face Fit Check: To check fit, place both hands along the edges of the respirator and exhale. If air escapes around your nose, adjust the nosepiece. If air leaks at the respirator edges, unfold or untwist the chin piece, repeat Steps 5 and 6.



8. PFL Mask correctly fit to face.

Before occupational use of this respirator, a written respiratory protection program must be implemented meeting all the local government requirements. In the United States, employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training and fit testing.

Warning: This product greatly reduces but does not eliminate the risk of transmission of infectious diseases. It is a single-use respirator; dispose of it properly after one use, and immediately if it is contaminated with blood or other bodily fluids.

AlphaProTech Not made with Natural Rubber Latex

Salt Lake City, Utah 84116

1.800.749.1363



Alpha Pro Tech, Inc.
 236 N. 2200 West
 Salt Lake City, Utah 84116
 1-800-749-1363



This respirator is approved only in the following configuration:

TC-	Protection ¹	Respirator	Cautions and Limitations ²
	N95	695	
84A-0457		X	A B C J M N O P

¹Protection:

N95 Particulate Filter (95% filter efficiency level) effective against particulate aerosols free of oil; time use restrictions may apply.

²Cautions and Limitations:

- A Not for use in atmospheres containing less than 19.5% oxygen.
- B Not for use in atmospheres immediately dangerous to life or health
- C Do not exceed maximum use concentrations established by regulatory standards.
- J Failure to properly use and maintain this product could result in injury or death.
- M All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA, and other applicable regulations.
- N Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- O Refer to user instructions and/or maintenance manuals for information on use and maintenance of these respirators.
- P NIOSH does not evaluate respirators for use as surgical masks.

NIOSH Approved

N-95 Particulate Respirator
Positive Facial Lock®
 With Magic Arch®



Cut-away View

Part No. 695
Contents: 35

*Fit Test required per NIOSH

- Patented Positive Facial Lock®
 "One Size Fits All"
- Unique elastic retention bands for easy donning
- Fluid-resistant
- Soft comfortable materials throughout

U.S. Patent #3985132, 4688566

AlphaProTech

PFL®
N-95 Particulate Respirator

• Positive Facial Lock® (PFL®)
 Allows total facial closure for "custom fit" on any size face when donned correctly.

• Integral Magic Arch®
 Creates a large comfortable breathing chamber with no restriction of downward vision.

• Traditional Pleated Design
 Saves inventory space, cannot be crushed when stored, and is easy to carry with you.

Test results for filtration available upon request.

GLUE SEAM

GLUE AREA

HIDDEN WHEN SET-UP

HIDDEN WHEN SET-UP

GLUE AREA