

## User Instruction

For special fitting demand, N95 respirator MS6815S, TC-84A-6810



### Warning

This respirator helps protect against certain particles. **Misuse may result in sickness or death.** For proper use, see supervisor, *User Instructions*, or call +86-512-66100068

Suzhou Sanical Protective Product Manufacturing Co., Ltd.  
No.6, Xihenggang Street, Yangchenghu Town, Xiangcheng District, Suzhou, China  
Tel: +86-512-66100068

### USE LIMITATIONS:

- Do **NOT** use the respirator or enter or stay in a contaminated area under the following circumstances:
  - *Atmosphere contains less than 19.5% oxygen*
  - *Atmosphere contains oil aerosols if using N protection class respirator*
  - *For protection against gases or vapors*
  - *Contaminants or their concentrations are unknown or immediately dangerous to life or health*
  - *Concentrations or contaminants exceed maximum use concentrations in applicable OSHA standards or applicable government regulations or 10 times the PEL (Permissible Exposure Limit), whichever is lower*
  - *For sandblasting, paint-spray operations, asbestos*
- Do **NOT** modify or misuse the respirator.
- Do **NOT** use the respirator with beards or other facial hair that interferes with direct contact between the face and the edge of the respirator, or any other conditions that may prevent a good face-seal.

### WARNING:

This respirator helps protect against certain contaminants. It may not eliminate the risk of contracting disease or infection. **Misuse may result in sickness or death.** Materials which may come into contact with the wearer's skin could cause allergic reactions to susceptible individuals. For more information or assistance please contact with Manufacturer.

### FITTING INSTRUCTIONS:



Follow these instructions each time respirator is worn.

- Hold the respirator in one hand, nose bridge at your fingertips, allowing the headbands to hang below your hand.
- Press the respirator against your face with the nose bridge of the respirator fit well with your nose
- Place the top band high on the back of your head. Move the bottom band over your head and position it below your ears.
- Fit-check the face-seal as follows:
  - *Place both hands over the respirator without disturbing its position.*
  - *Exhale sharply, if you detect air-leaks, if air leaks around your nose, replace the upper band; if air leaks around the other edges, reposition the lower band for better fit.*

**WARNING:** If you cannot achieve satisfactory fit, do **NOT** enter the contaminated area, you can consult manufacturer

### USE INSTRUCTIONS:

- Failure to follow all instructions and limitations concerning the use of the respirator and/or failure to wear the respirator for the whole duration of exposure to contaminants can seriously reduce the performance characteristics of the respirator and lead to illness, injury or death.
- Before use, the wearer must first be trained by the employer in the correct use of the respirator in accordance with applicable safety and health standards.
- OSHA standards including 29CFR1910.134 (e) (5) and CSA standard Z94.4-93 require that the wearer be fit-tested.
- Discard the respirator and replace with a new one if:
  - *Excessive clogging of the respirator causes breathing difficulty.*
  - *The respirator becomes damaged.*
- Leave the contaminated area if dizziness, irritation or other distress occurs.