

The NIOSH approval label will be printed on the box.

**JFY44351 NIOSH
P95 Particulate Respirator
With exhalation valve**

User's Instruction Manual

**JINFUYU INDUSTRIAL CO., LTD.
31F-9, NO. 97, Sec. 1, Xintai 5th Rd. Xizhi District,
New Taipei City, Taiwan**



**P95: 95% filter efficiency level,
Effective against particulate aerosol**

**ATTENTION: READ APPLICATIONS, LIMITATIONS
AND FITTING INSTRUCTIONS ON BOX BEFORE
RESPIRATOR IS WORN.**

CONTENTS: 5 RESPIRATORS

APPLICATIONS

USE AGAINST:

NIOSH approved for protection against all particulate, appropriate for prolonged use in both oil and non-oil containing environments.

DO NOT USE AGAINST:

1. Gases, vapors, including those present in paint spraying operations, asbestos or sandblasting.
2. Any contaminant whose concentration is unknown, or immediately dangerous to life.
3. Concentrations of Dust and Mists which exceed the maximum use concentration or 10 times contaminant's permissible exposure limit established by OSHA.

WARNING:

1. These respirators do not supply oxygen. Use only in well ventilated areas containing sufficient oxygen to support life.
2. Leave contaminated area immediately if the respirator becomes damaged and proper fit cannot be obtained or breathing becomes difficult. Dispose and change the respirator. If used in environments containing only oil aerosols, dispose of respirator after 40 hours of use or 30 days, whichever is first.

IMPORTANT LIMITATIONS:

1. All instructions and warnings have to be followed. Failure to do so will reduce the efficiency of respirator and cause sickness or death.
2. Training on the proper use of respirator must be provided by the employer in accordance with the safety and health standard.
3. Before occupational use of this respirator a written respiratory protection program must be implemented meeting all the local government requirements. In the United States employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training, and fit testing.
4. Do not alter or modify the respirator.
5. Not suitable for persons with beard or other facial hair which may reduce the efficiency of respirator.
6. Store the respirators in cool place. Avoid direct exposure of sunlight during storage as the filter media in this respirator can be degraded by extended exposure to direct sunlight.

Instructions for the 44351 P95 Maintenance free respirator

Fitting Instructions

The respirator user must be fit tested according to OSHA regulation 29 CFR 1910.134.

1. Hold the respirator in your hand with the nose area at your fingertips, allowing the head straps to hang below your hand.
2. Press the respirator against your face with the nose area on the bridge of your nose.
3. Place the top strap high on the back of your head. Move the bottom strap over your head and position it below your ears.
4. Adjust the strap tension by pulling the straps as shown in the Figure.
5. Cup both hands over the respirator and inhale vigorously. If no negative pressure is felt inside the respirator, adjust and tighten the nosepiece. If any leakage is detected, readjust the position of the respirator according to steps #2 to #5.

If you CANNOT achieve a proper fit, DO NOT enter the contaminated area.
See your supervisor for proper fitting respiratory protection.

