

USE INSTRUCTIONS

1. Failure to follow all instructions and limitations concerning the use of the respirator and/or failure to wear the respirator for the whole duration of exposure to contaminants can seriously reduce the performance characteristics of the respirator and lead to illness, injury or death.
2. Before Occupational Use of this respirator, a written respiratory protection program must be implemented meeting all the local government requirements. In United States, employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training, and fit testing.
3. Discard the respirator and replace with a new one if:
 - Excessive clogging of the respirator causes breathing difficulty.
 - The respirator becomes damaged.
4. Visually inspect the filtering facepiece prior to each use to assure that there are no holes in the breathing zone, especially around the weld attachment of the headbands to insure no damage has occurred. Discard respirator if any visual signs of holes or damage.
5. Leave the contaminated area if dizziness, irritation or other distress occurs.

USE LIMITATIONS

1. Do **NOT** use the respirator or enter or stay in a contaminated area under the following circumstances:
 - Atmosphere contains less than 19.5% oxygen.
 - Atmosphere contains oil aerosols.
 - For protection against gases or vapors.
 - Contaminants or their concentrations are unknown or immediately dangerous to life or health.
 - Concentrations or contaminants exceed maximum use concentrations in applicable OSHA standards or applicable government regulations or 10 times the PEL (Permissible Exposure Limit), whichever is lower.
 - For sandblasting, paint-spray operations, asbestos.
2. Do **NOT** modify or misuse the respirator.
3. Do **NOT** use the respirator with beards or other facial hair that interferes with direct contact between the face and the edge of the respirator, or any other conditions that may prevent a good face-seal.

WARNING

This respirator helps protect against certain contaminants. It may not eliminate the risk of contracting disease or infection. **Misuse may result in sickness or death.**

Materials which may come into contact with the wearer's skin could cause allergic reactions to susceptible individuals.

FOR MORE INFORMATION or assistance on Honeywell® products, contact your local Honeywell distributor.

Honeywell Mexicali

Laser 1199, Marán, 21385, Mexicali, B.C.

FITTING INSTRUCTIONS

Follow those instructions each time respirator is worn.

1. Cup the respirator in the hand, allowing the head straps to hang freely.
2. Hold the respirator under the chin with the nose-piece facing outwards. Place the lower head strap around the neck as seen in picture 2 above.
3. Holding the respirator against the face with one hand, place the top head strap above the ears, around the crown of the head.
4. Place both fingers on the nose clip and apply pressure onto the nose clip, forming it to the shape of your nose.
5. Place hands on each side of the respirator and move slightly right left, up and down to adjust the position of the respirator and achieve the most optimal fit on the face.

FIT-CHECK THE FACE-SEAL AS FOLLOWS

1. Place both hands over the respirator – do not disturb the position of the respirator;
2. Inhale sharply (a negative pressure should be felt inside the respirator);
3. If you detect air-leaks, readjust the head straps and/or the nose-piece.

WARNING:

If you cannot achieve satisfactory fit, do **NOT** enter the contaminated area. See your supervisor.

STORAGE

1. The respirator should be stored in an area which is well-ventilated, avoid direct sun exposure and dry area. Keep away from fire and pollution.
2. Storage temperature range is -20 °C (-4 °F) to 38 °C (100.4 °F). Maximum Relative Humidity <80%.



Storage Temperature
-20 °C (-4 °F) to
38 °C (100.4 °F)



Maximum
Relative
Humidity <80 %