

LUVMED® N95 Particulate Respirator

Model #: LUVMED 51-9501

Standard Facemask w/ Headband

CAUTION

All users must agree to fully review the following standard procedure guidelines BEFORE physically handling and/or wearing facemask respirators. Please retain this sheet for future reference. Additional important information from the supplier is included with each facemask.

BEFORE occupational use of this respirator, a written respiratory protection program must be implemented meeting all the local government requirements. In the United States, employers must comply with **OSHA 29 CFR 1910.134** which includes medical evaluation, training, and fit testing.

! IMPORTANT WARNING !

Improper use of this product can result in severe reactions such as sickness or death. Carefully follow all instructions as provided. Contact Info@AidwayCares.com with any questions regarding the correct use of masks.

PURPOSE / EFFECTIVENESS

The LUVMED N95 respirator mask is designed to filter **≥95%** particles that result from handling substances such as flour, pollen, wood, metals, minerals, silica, coal, ores, and others. Wear when performing such actions as sweeping, sanding, sawing, leaf blowing, welding, cutting, and metal heating. The masks can protect from droplets of non-oil aerosols or vapors emitted by sprays.

Do NOT Use to protect from oil aerosols, vapors or gases. Also not intended for use while sandblasting. The mask is ineffective against concentrated particles that are more than ten times the limit set by government regulations concerning occupational exposure.

REQUIRED SAFETY MEASURES

- Mask wear instructions must be understood and followed. A correct fit is mandatory.
- DO NOT Use mask in case of immediate danger where contaminants are already prevalent and concentrated.
- If the area is contaminated and breathing is inhibited or faintness occurs, withdraw from it immediately.
- If the mask suffers damage or it limits free breathing, dispose of it and use a new one.
- Before use, ensure the mask is free of flaws and is in good condition.
- Avoid touching, modifying, or washing mask over continued use.
- Not recommended for use with prevalent facial hair such as beards or mustaches that can keep the mask from adhering properly.
- If the mask incurs crumpled or creased damage, discard and replace.
- Use only per established regulation limits as set by your government.
- Use special caution during physical activity. Sudden contortion of the nose clip can lead to serious personal injury.
- All instructions must be followed. Heed all warnings. Wear only when particle exposure is likely, otherwise the mask's potency can cause illness or disability.

SPECIFICATIONS FOR STORAGE



LuvMed® N95 Respirators have a shelf life of three years from the manufacture date For maximum effectiveness they must be stored in an environment that is -4°F to 100°F (-20°C to +38°C) with 80% or less relative humidity).

DONNING AND DOFFING

1. Unfold the facemask with both hands, placing to your face, line up the metal nose clip horizontally at the bridge of your nose. Maintain that position and pull the mask's bottom over your chin.
2. Stretch the upper elastic strap over your head and position it a few inches above the ears.
3. Stretch the lower elastic strap over your head and position it below the ears, around the neck.
4. Using the index fingers of both hands, firmly press the pliable nose clip against the bridge of the nose and contour it by applying pressure while sliding the fingers down each side of the clip. Repeat to insure a tight fit against the nose.
5. Check to ensure the facemask is tightly sealed to the face by placing both hands directly over the mask, without affecting its positioning. Take a deep breath and exhale with force. There should be no air released from the bridge or mask edges.
6. In case of air leakage from the nose clip, repeat the molding process. If you are unable to achieve an tight seal, DO NOT USE and discard the mask.

REMOVING THE RESPIRATOR

1. DO NOT TOUCH the front of the respirator! It may be contaminated!
2. Remove by pulling the bottom strap over back of head, followed by the top strap, without touching the respirator.
3. Discard in waste container. Wash your hand!

	<p>Aidway Personal Care Product, Inc. 19465 E. Walnut Dr. N City of Industry, California, USA Toll Free: 1-866-624-3929</p>		
<p>This respirator is approved in only the following configuration</p>			
TC	Protection	Respirator	Cautions and Limitations
TC-84A-XXXX	N95	LUVMED 51-9501	ABCJMNOP
<p>1. Protection N95 - Particulate Filter (95% filter efficiency level). Effective against particulate aerosols free of oil; time use restrictions may apply.</p>			
<p>2. Cautions and Limitations A- Not for use in atmospheres containing less than 19.5% oxygen B- Not for use in atmospheres immediately dangerous to life or health C- Do not exceed maximum use concentrations established by regulatory standards. J- Failure to properly use and maintain this product could result in injury or death M- All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA, and other applicable regulations. N- Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer. O- Refer to User Instructions and/or maintenance manuals for information on use and maintenance of these respirators. P- NIOSH does not evaluate respirators for use as surgical masks.</p>			

Manufactured by: Aidway Personal Care Product Inc.
 19465 E Walnut Dr N., City of Industry, CA 91789 USA.
Toll free: 1-866-624-3929
LuvMed® is a trademark of Aidway Personal Care Product Inc.