



FIT-CHECK THE FACE-SEAL AS FOLLOWS:

- Place both hands over the respirator – do not disturb the position of the respirator.
- Inhale sharply (a negative pressure should be felt inside the respirator).
- If you detect air-leaks, readjust the head straps and/or the nose-piece.

The respirator must be used in accordance with OSHA regulations 29 CFR 1910.134 which includes the requirement for fit testing before use.

WARNING: If you cannot achieve satisfactory fit, do NOT enter the contaminated area. See your Supervisor or mail:yichita@yqgroup.vip to contact Shanghai Yunqing Industrial Co.,Ltd

Shanghai Yunqing Industrial Co.,Ltd..

2nd Buliding 1818 Shangs Road,Fengxian Distict,Shanghai.

Mail: yichita@yqgroup.vip

CAUTION: The straps on the respirator are completely latex-free. Individuals who are highly sensitive to natural rubber latex can be easy to use.

DONNING & FITTING INSTRUCTIONS:

Follow these instructions each time respirator is worn.

1. Cup respirator in one hand with nosepiece at fingertips, allow headbands to hang freely below hand.
2. Hold respirator under chin, with nosepiece up.
3. Locate the upper strap across the crown of the head and the lower strap below the ears. Straps must not be twisted.
4. Using both hands, mold noseclip to the shape of the nose to ensure a close fit and good seal. Pinching the noseclip using only one hand may result in less effective respirator performance.
4. Use two hands perform seal check prior to each wearing. To check fit, place both hands completely over the respirator and exhale. Be careful not to disturb the position of the respirator. If air leaks around nose, readjust the nosepiece as described in step. If air leaks at the respirator edges, work the straps back along the sides of your head. If you CANNOT achieve a proper fit, DO NOT enters the contaminated area. See your supervisor.