



1812

# Particulate Respirator N95

## User Instructions

### Important

Wearer must read fitting instructions and use limitations before use. Misuse of this respirator may result in sickness or death.

### Use For

Respirator with proper use protects against certain liquid or non-oil based particles from substances that do not emit oil vapors or aerosols. Respirator provides at least 95% filtration efficiency against solid and liquid aerosols that do not contain oil (see insert for approval label).

### Do Not Use For

Gases and vapors, oil aerosols, sandblasting, asbestos, lead, arsenic. The respirator does not supply oxygen.

### Service Time Limitation

All respirators should be replaced whenever they are damaged, soiled, or cause noticeably increased breathing resistance. Use and reuse of respirator subject only to considerations of hygiene, damage, and increased breathing resistance. Do not use the same respirator for more than 8 hours, continuous or intermittent, in dirty workplaces that could result in high filter loading. Leave contaminated area and remove and dispose of respirator if breathing becomes difficult or proper fit cannot be obtained.

### Warnings

This respirator does not eliminate the risk of contracting disease or infection. This product and its components are not made with natural rubber latex. Individuals with a compromised respiratory system should consult a physician and must complete a medical evaluation prior to use. Before occupational use of this respirator, a written respiratory protection program must be implemented meeting all the local government requirements. In the United States, employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training, and fit testing.

## ⚠ WARNING

User must follow all instructions and warnings on the proper use of this respirator. Misuse of this respirator may result in sickness and/or death. It is the user's responsibility to obtain a proper fit (see *Fitting Instructions*).

For additional information on this product please call 1-866-877-1381.

### United States Mask LLC.

4500 Westgrove Drive Suite 315  
Addison, Texas 75001

Proudly made in the USA



United States Mask LLC.  
4500 Westgrove Dr. Suite 315  
Addison, TX 75001  
1-866-877-1381



This respirator is approved only in the following configuration:

TC-	Protection	Respirator	Cautions and Limitations <sup>2</sup>
		1812	
84A-9393	N95	X	ABCJMNOP

### 1. Protection

N95 - Particulate Filter (95% filter efficiency level). Effective against particulate aerosols free of oil; time use restrictions may apply.

### 2. Cautions and Limitations

- A. Not for use in atmospheres containing less than 19.5% oxygen.
- B. Not for use in atmospheres immediately dangerous to life or health.
- C. Do not exceed maximum use concentrations established by regulatory standards.
- J. Failure to properly use and maintain this product could result in injury or death.
- M. All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA, and other applicable regulations.
- N. Never substitute modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- O. Refer to User Instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
- P. NIOSH does not evaluate respirators for use as surgical masks.



### FITTING INSTRUCTIONS

USER MUST FOLLOW THESE INSTRUCTIONS EACH TIME RESPIRATOR IS WORN.

1. Hold the bottom strap of the respirator using both hands, with the label facing away from you and the nosepiece facing up.
2. Place the respirator under your chin, stretch the bottom strap over your head, and place the strap below your ears.
3. Stretch and position the top band on the upper back of your head.
4. Using both hands, press the sides of the nosepiece to mold the wire to your nose.
5. Test the fit. To test the seal of your respirator 1) cup both hands over the respirator being careful not to disturb the position 2) Take several large exhales. If air leaks around the edges, reposition the straps or adjust strap tension for a better fit. **If you cannot achieve a proper seal, do not enter the contaminated area. See your supervisor.**