



Surgical N95 Respirator with Comfort Strip

Wearing Instructions

- Follow the **Directions for Application** (*below*) each time a respirator is worn.
- Perform a **Fit Check** prior to each wearing. To check the respirator-to-face seal, cover the front of the respirator completely with both hands, being careful not to disturb the position of the respirator. Inhale Sharply. A negative pressure should be felt inside the respirator. If any leakage is detected, readjust the position of the respirator according to step #8.
- If proper fit cannot be obtained, **DO NOT** enter a contaminated area.
- Wear during all times of exposure.

Directions for Application



1

Insert fingers to form respirator into a cup shape.



2

Hold the respirator upside-down to reveal the two headbands.



3

Separate the headbands with your index fingers and thumbs.



4

Cup the respirator under your chin while continuing to hold the headbands with your index fingers and thumbs.



5

Pull the headbands over your head.



6

Position the lower headband over the base of your head and



7

Position the top headband on the crown of your head.



8

Form the nosepiece to the bridge of your nose and adjust respirator until you have a secure fit.

Failure to follow all instructions and warnings on the use of this respirator, especially instructions for a proper fit, will reduce respirator effectiveness which may result in sickness or death.

DESCRIPTION: The Surgical/N95 Respirator is NIOSH approved as an N95 particulate filter respirator. It provides 95% filtration efficiency of 0.3 micron particles. It meets the CDC Guidelines for TB exposure control. **INTENDED USE:** The Surgical/N95 Respirator is intended for use by operating room personnel and other health care workers to protect both patients and health care workers from transfer of microorganisms and airborne particulate materials and reduce the potential exposure of the wearer to blood and body fluids. **CONTRAINDICATION:** This device is not designed, sold, or intended for use except as indicated. **SPECIAL OR CRITICAL USER INSTRUCTIONS:** 1. This respirator has been approved as a NIOSH N95 filtering facepiece respirator, for use in healthcare settings, as a Surgical N95 Respirator conforming to recognized standards for biocompatibility, flammability, and fluid resistance. 2. Before occupational use of this respirator, a written respiratory protection program must be implemented meeting all the local government requirements. In the United States, employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training, and fit testing. 3. If the respirator becomes damaged or breathing becomes difficult, leave the area. Discard and replace the respirator. **WARNING:** This respirator helps protect against certain particulate contaminants but does not eliminate exposure to or the risk of contracting any disease or infection. **REMINDER:** It is very important to fit check the respirator each time you don the respirator. Consult the facility's Respiratory Protection Program for further instruction.



Helpful Videos



Instrucciones de Uso
en Español



Instructions d'Utilisation
en Français

Made in U.S.A.

ASTM Level 3 (160 mm Hg) Fluid Resistance Protection

Manufactured by Prestige Ameritech
7201 Iron Horse Blvd
North Richland Hills, TX 76180
www.progearmasks.com
1-866-401-8972

NIOSH-029-03