



Fido disposable particulate respirator

**F720 N95 Particulate Respirator**

**User's Instructions**

**WARNING TO USER**

Refer to the user manual, technical and other all user instructions, warnings and the limitations of the respirator. This respirator provides limited protection and user is responsible and must be trained to a competent person to provide protection from various concentrations according to OSHA 29 CFR 1910.134 an applicable safety and health standard. The respirator's life and use must be checked each time before entering a contaminated area with the respirator used to ensure that all parts of the respirator are functioning. Failure to follow all of these instructions, warnings and limitations, contained in this label may cause the effectiveness of this respirator and result in serious sickness or DEATH.

The F720 N95 respirator with a filtration efficiency of at least 95% against solid or liquid particles that do not settle out.

**DO NOT USE FOR:**

The F720 N95 respirator should not be used for gases, vapors, oil mists or where concentrations of contaminants exceed OSHA or applicable safety and health standard's maximum use concentrations, or where a manufacturer's warnings are left and health.

**USE INSTRUCTIONS:**

- Before occupational use of this respirator a written respiratory protection program must be implemented meeting all local government requirements. In the United States, employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training, and fit testing.
- Before use, carefully inspect the entire respirator for correct assembly and for tears or damage including the straps, clips and adjustment parts (if applicable). **Discard the entire respirator if any part does not function properly and get a new respirator.**
- If the respirator is to be used in a contaminated area, immediately and per a new respirator in face and cannot be re-used, leave the contaminated area immediately and get a new respirator.
- If leakage, irritation, dizziness or other distress occurs, leave the contaminated area immediately and replace your respirator or another professional.
- Do not alter, abuse or misuse this respirator.
- Periodically inspect and maintain in accordance with applicable regulations.

**USE LIMITATIONS:**

- This respirator does not supply oxygen and must not be used in atmospheres containing less than 19.5% oxygen.
- This respirator is not to be used when the eyes are irritated, facial hair, scars or other conditions that prevent a good seal between the face and the respirator.
- This respirator is designed for occupational or professional use by adults who are properly trained in their use and limitations. Respirators may not be used by children.
- Particulates and other contaminants which may be dangerous to your health include those that you cannot see, taste or smell. Respirators can help to protect you from respirator cannot address contaminants that cannot be seen, tasted or smelled through other means such as the skin. Additional types of personal protective equipment may be required.

**STORAGE:**

- Store unused respirators in their sealed, original packaging in a clean, dry, non-contaminated area that is protected from extreme temperatures, moisture, sunlight, dust and chemicals.
- Never attempt to clean or reuse the respirator.
- Used respirators should not be stored and should be disposed of properly.



Fido disposable particulate respirator



N95



FILTERING FACEPIECE APPROVAL LABEL



FIDO MASKS CO., LTD.  
No. 8-5 She-Ling Li, Yuan-Li Town,  
Miao-Li County 35872, Taiwan  
TEL: 886-37-741471

**THIS RESPIRATOR IS APPROVED ONLY IN THE FOLLOWING CONFIGURATION:**

IC	Protection <sup>1</sup>	Respirator	Cautions and Limitations <sup>2</sup>
84A-7992	N95	F720	X ABCMNOP

**1. Protection**

N95 - Particulate Filter (95% filter efficiency level)  
Effective against particulate aerosols free of oil.  
Time use restrictions may apply.

**2. Cautions and Limitations**

- A - Not for use in atmospheres containing less than 19.5% oxygen.
- B - Not for use in atmospheres immediately dangerous to life or health.
- C - Do not exceed maximum use concentrations established by regulatory standards.
- J - Failure to properly use and maintain this product could result in injury or death.
- M - All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA and other applicable regulations.
- N - Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- O - Refer to User's Instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
- P - NIOSH does not evaluate respirators for use as surgical masks.



**FITTING INSTRUCTIONS:**

Before fitting the respirator, carefully inspect the entire respirator to ensure that the respirator is not damaged and is functioning correctly. **Discard the entire respirator if any part does not function properly.**

- Hold the respirator in your hand with the nosepiece at your fingertips and allowing the straps to fall below the hand.
- Place respirator under the chin with the nosepiece (narrow end) up. Raise the continuous top strap over your head and place it high to the back of the head.
- Pull the continuous bottom strap over your head and position it around the neck below the ears.
- Using both hands, push the nosepiece inwards with your fingertips to mold and shape the nose area to the face.
- A seal check MUST be performed each time the user enters a contaminated area. Lightly cup both hands over the front of the respirator and INHALE SHARPLY. A negative pressure should be felt inside of the respirator. If any air flow is detected around the edges of the respirator, readjust the straps. Repeat seal check until a proper seal is achieved; otherwise, see your supervisor. **NEVER enter a contaminated area without a proper seal.**

**Note:** Pinching the nosepiece using one hand may result in improper fit and less effective respirator performance. Use two hands.

**Manufacturer Address:**

Fido Mask Co., Ltd  
8-5 She-Ling, Yuan-Li Town, Miao-Li County 35872, Taiwan

Produce date: YYYY/MM

Lot No. XXXXX

(the Lot No. XXXXX will be as same as abbreviated Label indicated, the "XXXXX" should be printed as purchase order number )

Contents: 20 pcs.

Fido disposable particulate respirator



N95