



S6N



Fido disposable particulate respirator

F720C N95 Particulate Respirator

User's Instructions

WARNING TO USER

Before use, the user must read additional and other all user instructions, warnings and the function of this respirator. The user must understand the level of protection and must be properly selected and fit tested to a competent person to avoid overexposure to hazardous contaminants as well as OSHA 29 CFR 1910.134 or applicable safety of health method. The respirator face seal must be checked on the user face to ensure a complete seal. The respirator must be worn **strictly** at all times of exposure to contaminants. Failure to follow all of these instructions, warnings and limitations, contained on this tag may reduce the effectiveness of this respirator and result in a serious risk to the user's health.

The F720C N95 respirator has a filtration efficiency of at least 95% against solid or liquid particles that do not settle out.

DO NOT USE FOR:

The F720C N95 respirator should not be used for gases, vapors, oil aerosols or other concentrations of contaminants except OSHA or applicable safety and health standard. Contaminants are concentrations, are unknown or immediately dangerous to life and health.

USE INSTRUCTIONS:

- Before respiratory use of this respirator in a certain environment protection program must be implemented according to the local government arrangements. In the United States, employers must comply with OSHA's 29 CFR 1910.134 which includes additional selection, testing, and fit testing.
- Proper use includes: carefully inspect the respirator for correct assembly and the type or change according to steps. This will eliminate valve of leakage with use. **Discard the entire respirator if any part does not function properly or use one respirator.**
- If the respirator becomes damaged, distorted, soiled, or breathing becomes difficult or a proper seal cannot be maintained, leave the contaminated area immediately and use a new respirator.
- If electronic indicators, sensors or other devices occur, leave the contaminated area immediately and replace your respirator or modify professional.
- Do not alter (trim or shape) the respirator.
- Properly dispose of used respirator in accordance with applicable regulations.

USE LIMITATIONS:

- This respirator does not supply oxygen and must not be used in atmospheres containing less than 19.5% oxygen.
- This respirator is not to be used when the user has head, facial hair, sores or other conditions that prevent a good seal between the face and the respirator.
- This respirator is designed for occupational or professional use by adults who are properly trained in their use and maintenance. Respirators must not be used by children.
- Facemasks and other contaminants which can be dangerous to your health include those that you cannot see, taste or smell. Respirators are held in place by a strap system which includes adjustable components. Use cannot prevent entry through other routes such as the skin. Additional types of personal protective equipment may be required.

STORAGE:

- Store unused respirators in their sealed, original packaging in a clean, dry, non-contaminated area that is protected from extreme temperatures, excessive sunlight, dust and dirt.
- Never attempt to clean and reuse this respirator.
- Used respirators should not be stored and should be disposed of properly.



See information supplied by the manufacturer

End of Strap Use 5 hours

Temperature range of storage conditions -10°C to 40°C

Maximum relative humidity of storage condition +75%

Fido disposable particulate respirator



FILTERING FACEPIECE APPROVAL LABEL



FIDO MASKS CO., LTD.
No. 8-5 Shee-Ling Li, Yuan-Li Town,
Miao-Li County 35872, Taiwan
TEL: 886-37-741471

THIS RESPIRATOR IS APPROVED ONLY IN THE FOLLOWING CONFIGURATION:

TC*	Protection ¹	Respirator	Caution and Limitations ²
84A-7793	N95	X	AHC/MNOP

I. Protection

N95 - Particulate Filter (95% filter efficiency level)
Effective against particulate aerosols free of oil;
tube use restrictions may apply

2. Caution and Limitations

- A - Not for use in atmospheres containing less than 19.5% oxygen.
- B - Not for use in atmospheres immediately dangerous to life or health.
- C - Do not exceed maximum use concentrations established by regulatory standards.
- J - Failure to properly use and maintain this product could result in injury or death.
- M - All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA and other applicable regulations.
- N - Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- O - Refer to User's Instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
- P - NIOSH does not evaluate respirators for use as surgical masks.



FITTING INSTRUCTIONS:

Before fitting the respirator, **carefully inspect the entire respirator** to ensure that the respirator is not damaged and is functioning correctly. **Discard the entire respirator if any part does not function properly.**

- Hold the respirator in your hand with the nosepiece at your fingertips and allowing the straps to fall below the hand.
- Place respirator under the chin with the nosepiece (narrow end) up. Pull the continuous bottom strap over your head and position it around the neck below the ears.
- Raise the continuous top strap over your head and place it high to the back of the head.
- Using both hands, push the nosepiece inwards with your fingertips to mold and shape the nose area to the face.
- A seal check **MUST** be performed each time the user enters a contaminated area. Lightly cup both hands over the front of the respirator and **INHALE SHARPLY**. A negative pressure should be felt inside of the respirator. If any air flow is detected around the edges of the respirator, readjust the straps. Repeat seal check until a proper seal is achieved; otherwise, see your supervisor. **NEVER enter a contaminated area without a proper seal.**

Note: Pinching the nosepiece using one hand may result in improper fit and less effective respirator performance. Use two hands.

Manufacturer Address:

Fido Mask Co., Ltd
8-5 Shee-Ling, Yuan-Li Town, Miao-Li County 35872, Taiwan

Produce date: YYYY/MM

Lot No. XXXXXX

(the Lot No. XXXXXX will be as same as abbreviated Label indicated, the "XXXXX" should be printed as purchase order number)

Contents: 20 pcs.

Fido disposable particulate respirator



N95

