



## Ford N95 Respirator

Ford Motor Company  
One American Rd.  
Dearborn, MI 48126  
Email: ppe1@ford.com



This respirator is approved only in the following configuration:

TC-	Protection <sup>1</sup>	Respirator	Cautions and Limitations <sup>2</sup>
		FMCN95	
84A- PHXX	N95	X	ABCJMNO PH

### 1. Protection

N95 - Particulate Filter (95% filter efficiency level). Effective against particulate aerosols free of oil; time use restrictions may apply.
---

### 2. Cautions and Limitations

- A Not for use in atmospheres containing less than 19.5% oxygen.
- B Not for use in atmospheres immediately dangerous to life or health.
- C Do not exceed maximum use concentrations established by regulatory standards.
- J Failure to properly use and maintain this product could result in injury or death.
- M All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA, and other applicable regulations.
- N Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- O Refer to User Instructions and/or maintenance manuals for information on use and maintenance of these respirators.
- P NIOSH does not evaluate respirators for use as surgical masks.
- PH The limited FFR approval will be effective until the declaration that circumstances exist justifying the authorization of the emergency use of personal respiratory protective devices during the COVID-19 outbreak is terminated under Section 564(b)(2) of the Act, or the EUA is revoked under Section 564(g) of the Act, or until such time that the approved configuration can no longer be maintained.

Ford N95 Respirator (PHF0615)

# N95 Respirator User Instructions

## Installation Process

- Step 1:** Remove the respirator from its packaging and hold with straps facing upward.
- Step 2:** Fully open the top and bottom panels, holding at the center of the nose cushion. Verify that the top and bottom panels are unfolded and completely opened.
- Step 3:** Place the respirator on your face so that the foam nose cushion rests on your nose and the bottom panel is securely under your chin.
- Step 4:** Pull the top strap over your head, positioning it high on the back of the head. Next, pull the bottom strap over the head and position it around the neck and below your ears.
- Step 5:** Adjust the respirator for a comfort by pulling the top panel toward the bridge of your nose and the bottom panel under your chin. Verify that hair, facial hair, jewelry and/or clothing are not between your face and the respirator as they will interfere with fit.
- Step 6:** Place fingertips from both hands at the top of the metal nose arch band. Using two hands, press the index finger on the nose arch band to mold it to the shape of the nose, pushing inward while moving your fingertips down both sides of the nose arch band. Always use two hands when molding the nose arch band. Do not pinch the nose arch band with one hand, as this may result in improper fit and less effective respirator performance.
- Step 7:** Verify the seal of the mask to the face. Place one or both hands completely over the middle panel of the respirator. Inhale and exhale sharply. Be careful not to disturb the position of the respirator. If air leaks around your nose, re-adjust the nose arch band as described in Step 6. If air leaks around respirator edges, adjust panels and the position of the straps, verifying the respirator edges fit snugly against the face. If you cannot achieve a proper seal, do not enter a contaminated area.



**Before occupational use of this respirator, a written respiratory protection program must be implemented meeting all the requirements of OSHA 29 CFR 1910.134 such as training, fit testing, medical evaluation, and applicable OSHA substance specific standards.**

## N95 Respirator User Instructions

### Removal Process

*Note: The removal process can be completed with one or two hands.*

**Step 1:** Without touching the respirator facepiece, slowly lift the bottom strap from around your neck up over your head.

**Step 2:** Lift off the top strap. Do not touch the respirator.

**Step 3:** Disposition the respirator according to the infection control policy.



**Before occupational use of this respirator, a written respiratory protection program must be implemented meeting all the requirements of OSHA 29 CFR 1910.134 such as training, fit testing, medical evaluation, and applicable OSHA substance specific standards.**

Ford N95 Respirator (PHF0615)