1740c1

550489 rev L

550489L-1740c1DISP-030609mDR3.pdf

5.25"

1740 N95 Respirator
With Exhalation Valve

1740 N95 Respirator
With Exhalation Valve

LO

84A-160	N95	1730	1740	
84A-161	N95		X	ABCJMNOP

- Protection
 - N-95 Particulate filter (95% filter efficiency level) effective against particulate aerosols free of oil; time use restrictions may apply.
- Cautions and Limitations
 - A - Not for use in atmospheres containing less than 19.5% oxygen.
 - B - Not for use in atmospheres immediately dangerous to life or health.
 - C - Do not exceed maximum use concentrations established by regulatory standards.
 - J - Failure to properly use and maintain this product could result in injury or death.
 - M - All approved respirators shall be selected, fitted, used and maintained in accordance with MSHA, OSHA and other applicable regulations.
 - N - Never substitute, modify, add or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
 - O - Refer to user instructions and/or maintenance manuals for information on use and maintenance of the respirators.
 - P - NIOSH does not evaluate respirators for use as surgical masks.

IMPORTANT NOTICE TO PURCHASER
 Manufacturer warrants this product to be free from defects in material and workmanship. In the event that this product proves upon inspection to contain defects in material or workmanship, Manufacturer undertakes, at its option, to replace any portion proving defective or to refund to the Buyer the purchase price paid by Buyer for this product. This warranty and Buyer's remedy of replacement or refund constitute the Manufacturer's exclusive warranty and the Buyer's sole remedy, in lieu of all other warranties expressed or implied. Except as provided above, Manufacturer shall not be liable for any loss or damage, whether direct, consequential, or incidental occasioned by the sale, use or misuse of this product, other than death or personal injury directly resulting from Manufacturer's negligence. Buyer shall determine the suitability of this product for its intended use and except provided above, assumes all risks and liability in connection therewith.

LOUIS M. GERSON CO., INC.
 Middleboro, MA 02346 U.S.A.
 © 2009 All Rights Reserved. Made in U.S.A.
 550489 rev L
 www.GersonCo.com

WARNINGS

Read and understand before use.

- A properly selected respirator is essential to protect your health. Before using this or any product consult an Industrial Hygienist or Occupational Safety Professional to determine its suitability for your intended use.
- This product will not protect against gases or vapors. Do not use in spray paint operations.
- Do not use in or around sand blasting.
- Do not use for protection against asbestos.
- This product does not supply oxygen. Use this respirator in a well ventilated area containing sufficient oxygen to support life.
- Do not use when concentrations of contaminants are unknown or immediately dangerous to life.
- Leave work area immediately if:
 - a) breathing becomes difficult, or
 - b) dizziness or other signs of distress occur.
- If user has any illness or disorders, consult a licensed medical doctor to determine suitability of use.
- Facial hair and certain facial characteristics may prohibit effective use of this product.
- Use only in accordance with instructions, labels and limitations relating to this product.
- Do not alter or modify this product in any way.
- Follow the manufacturer's instructions for fitting the respirator to the face.
- Filtering facepieces are to be inspected prior to each use to assure there are no holes in the breathing zone and that no damage has occurred.

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Cat. No. 081740

1740 N95 Respirator

NIOSH Approved N95 Particulate Respirator with Exhalation Valve

A Cooler More Comfortable Experience

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NOTE: Before issuing any respirator to be worn, a face fit test must be performed. See OSHA regulations 1910.134.

6.50"