

USER´S INSTRUCTIONS LIFE 1095 / LIFE1095-2
RESPIRATORS
(steps 1 - 2)

LABEL FOR BOXES ([step 1](#))



Calle 107 # 37-39
PBX: + 57 5 382 16 16
e-mail : ingenieriatecnica@indusaver.com
Barranquilla Colombia
NIOSH Approved – See the Bags

Read the instructions in the bags to
individual respirator

LABEL FOR INDIVIDUAL BAGS (step 2)

ADJUSTMENT INSTRUCTIONS

**Fitting instructions to be followed each time
the the respirator is:**



1. Hold the respirator in your hand, allowing the headbands to hang freely below your hands.



2. Position the respirator comfortably under your chin and against the bridge of your nose.



3. Pull the top headbands over your head placing it at the crown of your head. Pull the bottom headband over your head and place it below your ears.



4. USING BOTH HANDS, mold the metal nosepiece to the shape of your nose and cheek area. Using one hand to pinch the nosepiece may cause a bad fit, resulting in lower effective respirator performance.

5. FIT CHECK: Place both hands completely over the respirator. Inhale sharply. Be careful not to disturb the position of the respirator. A negative pressure (mask tends to collapse) should be felt. If not, or if air leaks around the edges of the respirator and headbands for a better fit. If unable to adequately seal filter with hands, obtain the optional respirator fit check cup to perform this test. If you CANNOT achieve a proper fit, DO NOT WEAR THIS MASK.



NIOSH APPROVED

LIFE 1095 - N-95

TC 84A-5350

CAUTION



This respirator helps to protect against certain particulates. The improper use can cause illness or death.

“Before occupational use of this respirator, a written respiratory protection program must be implemented meeting all the local government requirements. In the United States, employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training, and fit testing”