

JFY 668150, 668151 (N95)

Donning Procedures

1. Hold the respirator in one hand with the nosepiece at your fingertips, pulling the headbands by the other hand.
2. Press the respirator against your face with the nosepiece on the bridge of your nose.
3. Place the top band high on the back of your head. Move the bottom band over your head and position it below your ears.
4. Adjust the strap tension by pulling the ring adjuster as shown in the Figure.
5. Using both hands, mould the nosepiece to the shape of your nose.
6. Test the fit. Cup both hands over the respirator and exhale vigorously. If air flows around your nose, tighten the nosepiece. If air leaks around the edges, reposition the bands for better fit.

If you CANNOT achieve a proper fit, DO NOT enter the contaminated area.
See your supervisor for proper fitting respiratory protection.

