

Gateway Safety 80602V (R95)

Donning Procedures

1. Hold mask in one hand with straps around back of hand. Place over mouth and nose with nose clip over bridge of nose. Pull lower strap over head and fit around neck. Pull upper strap over head and place it above ears.
2. Form nose clip around nose and adjust tension in straps by positioning straps on head until you obtain a positive face seal.
3. Perform a User Seal Check prior to each wearing. It is the wearer's responsibility to obtain a proper fit. Each time respirator is used, wearer must test fit by cupping both hands over the respirator and inhaling vigorously. If no negative pressure is felt inside the respirator, adjust and tighten the nosepiece and reposition straps. If any leakage is detected, readjust the position of the respirator according to steps 2 and 3 until you get a better face seal. If you CANNOT achieve a proper face seal, DO NOT enter the contaminated area. See your supervisor.

