

JFY 44351 (P95)

Donning Procedures

1. Hold the respirator in your hand with the nose area at your fingertips, allowing the head straps to hang below your hand.
2. Press the respirator against your face with the nose area on the bridge of your nose.
3. Place the top strap high on the back of your head. Move the bottom strap over your head and position it below your ears.
4. Adjust the strap tension by pulling the straps as shown in the Figure.
5. Cup both hands over the respirator and inhale vigorously. If no negative pressure is felt inside the respirator, adjust and tighten the nosepiece. If any leakage is detected, readjust the position of the respirator according to steps #2 to #5.

If you CANNOT achieve a proper fit, DO NOT enter the contaminated area.
See your supervisor for proper fitting respiratory protection.

