



Manufactured by
PandMedic Inc.
4544 W. Russel Rd Ste.D
Las Vegas NV 89118



Contact:
1(800) 624-3129

This respirator is approved only in the following configuration:

TC-	Protection	Respirator	Cautions and Limitations
		MedicPro N95	
84A-PHnm	N95	X	ABCJMNO PH

1. Protection

N95 - Particulate Filter (95% filter efficiency level).
Effective against particulate aerosols free of oil;
time use restrictions may apply.

2. Cautions and Limitations

- A-** Not for use in atmospheres containing less than 19.5% oxygen.
- B-** Not for use in atmospheres immediately dangerous to life or health.
- C-** Do not exceed maximum use concentrations established by regulatory standards.
- J-** Failure to properly use and maintain this product could result in injury or death.
- M-** All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA, and other applicable regulations.
- N-** Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- O-** Refer to User Instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
- P-** NIOSH does not evaluate respirators for use as surgical masks.
- PH-** The limited FFR approval will be effective until the declaration that circumstances exist justifying the authorization of the emergency use of personal respiratory protective devices during the COVID-19 outbreak is terminated under Section 564(b)(2) of the Act, or the EUA is revoked under Section 564(g) of the Act, or until such time that the approved configuration can no longer be maintained.



MedicPro N95



MEDICAL GRADE N95
N95 PARTICULATE RESPIRATOR

PROUDLY MADE IN USA | 1 PC

Instructions



1
Unfold the mask and pull the bands outward on both ends.



2
Position the mask on your face with adjustable metal strip facing up.



3
Place the bands over the head



4
Press down on the metal strip lightly to ensure the mask is securely sealed.

Proper Care

Before occupational use of this respirator, a written respiratory protection program must be implemented meeting all the local government requirements. In the United States, employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training, and fit testing, and applicable OSHA substance-specific standards.

Warning: DO NOT use in areas with oxygen concentration below **19.5%**

This Mask Has N95 Filtration Capabilities

BF95

MADE IN USA



0 607462 089552

