



1836

Particulate Respirator N95

User Instructions

Important

Wearer must read fitting instructions and use limitations before use. Misuse of this respirator may result in sickness or death.

Use For

Respirator with proper use protects against certain liquid or non-oil based particles from substances that do not emit oil vapors or aerosols. Respirator provides at least 95% filtration efficiency against solid and liquid aerosols that do not contain oil (see insert for approval label).

Do Not Use For

Gases and vapors, oil aerosols, sandblasting, asbestos, lead, arsenic. The respirator does not supply oxygen.

Service Time Limitation

All respirators should be replaced whenever they are damaged, soiled, or cause noticeably increased breathing resistance. Use and reuse of respirator subject only to considerations of hygiene, damage, and increased breathing resistance. Do not use the same respirator for more than 8 hours, continuous or intermittent, in dirty workplaces that could result in high filter loading. Leave contaminated area and remove and dispose of respirator if breathing becomes difficult or proper fit cannot be obtained.

Warnings

This respirator does not protect against the risk of contracting disease or infection. Individuals who are highly sensitive to natural rubber latex may have an allergic reaction. Individuals with a compromised respiratory system should consult a physician and must complete a medical evaluation prior to use. The wearer should be trained in the proper use, including fit testing, in accordance with applicable safety and health standards per OSHA guideline 29 CFR 1910.134. If you cannot achieve a proper seal, DO NOT enter the contaminated area.

WARNING

User must follow all instructions and warnings on the proper use of this respirator. Misuse of this respirator may result in sickness and/or death. It is the user's responsibility to obtain a proper fit (see *Fitting Instructions*).

For additional information on this product please call 1-866-877-1381.

United States Mask LLC.

4500 Westgrove Drive Suite 315
Addison, Texas 75001

Proudly made in the USA



Approval: TC-84A-PHXX
Assigned to: United States Mask LLC.
4500 Westgrove Dr. Suite 315
Addison, TX 75001
1-866-877-1381



Filtering Facepiece Respirator for Public Health Emergency
Protection:

Read and understand UI prior to use.
User must adhere to all C&L in the UI.

The following Cautions and Limitations spelled out in the UI apply:
A, B, C, J, L, M, N, O, P, PH

CAUTIONS AND LIMITATIONS

- A.** Not for use in atmospheres containing less than 19.5 percent oxygen.
- B.** Not for use in atmospheres immediately dangerous to life or health.
- C.** Do not exceed maximum use concentrations established by regulatory standards.
- J.** Failure to properly use and maintain this product could result in injury or death.
- M.** All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA, and other applicable regulations.
- L.** Follow the manufacturer's User Instructions for changing cartridges, canister and/or filters
- N.** Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- O.** Refer to User Instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
- P.** NIOSH does not evaluate respirators for use as surgical mask.
- PH.** The limited FFR approval will be effective until the declaration that circumstances exist justifying the authorization of the emergency use of personal respiratory protective devices during the COVID-19 outbreak is terminated under Section 564(b)(2) of the Act, or the EUA is revoked under Section 564(g) of the Act, or until such time that the approved configuration can no longer be maintained.



FITTING INSTRUCTIONS

USER MUST FOLLOW THESE INSTRUCTIONS EACH TIME RESPIRATOR IS WORN.

1. Hold the bottom strap of the respirator using both hands, with the label facing away from you and the nosepiece facing up.
2. Place the respirator under your chin, stretch the bottom strap over your head, and place the strap below your ears.
3. Stretch and position the top band on the upper back of your head.
4. Using both hands, press the sides of the nosepiece to mold the wire to your nose.
5. Test the fit. To test the seal of your respirator 1) cup both hands over the respirator being careful not to disturb the position 2) Take several large exhales. If air leaks around the edges, reposition the straps or adjust strap tension for a better fit. **If you cannot achieve a proper seal, do not enter the contaminated area. See your supervisor.**