

Suzhou Sanical Protective Product Manufacturing Co., Ltd.

Respirator P385 (P95)

Donning Procedures

1. Open the respirator and hold in one hand with the nosepiece at your fingertips, allowing the headbands to hang below your hand.
2. Press the respirator against your face with the nosepiece on the bridge of your nose. Place the top band high on the back of your head. Move the bottom band over your head and position it below your ears.
3. The length of the headband can be adjusted. In order to better protect the protective effect, hold both the left and right ends of the upper or lower headbands with both hands and tighten them to the sides of the ears so that the mask fits perfectly on the face.
4. Place your two fingertips from both hands at the top of the nosepiece. Using two hands mold the nose area to the shape of your nose by pushing inward while moving your fingertips down both sides of the nosepiece.

Warning: Pinching the nosepiece using one hand may result in improper fit and less effective respirator performance. Use two hands.

5. Fit-check the face-seal as follows:

- Place both hands over the respirator without disturbing its position;
- Inhale sharply, a negative pressure should be felt inside the respirator;
- If a negative pressure doesn't be felt, there is air leak, please re-adjust the nosepiece fellow step 4, and/or, re-adjust head straps length fellow step 3.

Note: Do not wear with beards or other facial hair that may inhibit contact between the face and the edge of the respirator;

WARNING: If you cannot achieve satisfactory fit, do NOT enter the contaminated area. See your Supervisor.

