



Master & Frank Enterprise Co., Ltd.

19F, No. 57, Sec. 2, Tun Hwa S. Rd., Taipei City, Taiwan

Dongguan Shin Yi Healthcare Products Factory (Factory)

No. 17 Ban Hu Road, Huang Jiang Town, Dong Guan City, Guang Dong Province, P.R. China

86-769-83361381



THIS RESPIRATOR IS APPROVED ONLY IN THE FOLLOWING CONFIGURATION:

| TC- | Protection ¹ | Respirator | Cautions and Limitations ² |
|----------|-------------------------|------------|---------------------------------------|
| | | MF-0106 | |
| 84A-9527 | Surgical/N95 | X | ABCJM NOS |

1. Protection

N95 - Particulate Filter (95% filter efficiency level)

Effective against particulate aerosols free of oil;

time use restrictions may apply

2. Cautions and Limitations

A - Not for use in atmospheres containing less than 19.5% oxygen.

B - Not for use in atmospheres immediately dangerous to life or health.

C - Do not exceed maximum use concentrations established by regulatory standards.

J - Failure to properly use and maintain this product could result in injury or death.

M - All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA and other applicable regulations.

N - Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.

O - Refer to user instructions, and/or maintenance manuals for information on use and maintenance of these respirators.

S - Special or critical User Instructions and/or specific use limitations apply. Refer to User Instructions before donning.

S – Special or Critical User Instructions

This respirator has been approved as a NIOSH N95 filtering facepiece respirator, for use in healthcare settings, as a Surgical N95 Respirator conforming to recognized standards for biocompatibility, flammability, and fluid resistance at 160mmHg.

Note :

- 1 - Not for use in atmospheres containing less than 19.5% oxygen.
- 2 - Not for use in atmospheres immediately dangerous to life or health.
- 3 - Do not exceed maximum use concentrations established by regulatory standards.
- 4- Failure to properly use and maintain this product could result in injury or death.
- 5 - All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA and other applicable regulations.
- 6 - Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- 7 - Refer to user instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
- 8 - Special or critical User Instructions and/or specific use limitations apply. Refer to User Instructions before donning.
- 9- Before occupational use of this respirator a written respiratory protection program must be implemented meeting all the local government requirements. In the United States employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training, and fit testing.
- 10- Storage conditions: place in a dry and well ventilated place. Do not place directly in the sun.
- 11- Please pay attention to the expiration date on the bags or box before use. Use of the product, over the expiration date is strictly prohibited.
- 12- The respirators are disposable. The respirators cannot be cleaned or reused.
- 13- Do not use near fire, flames or extremely hot objects.
- 14- Please check the respirator surface for damage. Make sure the nosepiece and headstraps are firmly attached before opening the package. Do not use if there is damage.
- 15- Please refer to the user instructions for use when wearing.
- 16- If you cannot achieve a proper fit, DO NOT enter the contaminated area. See your supervisor or call 0086-0769-83361381 Dongguan Shin Yi Healthcare Products Factory.

Donning & Fitting Instruction



1. Hold respirator in position over nose and mouth, with retaining straps around back of hand.



2. Pull the lower strap over head and fit around neck. Pull the upper strap over head and fit above ears.



3. Adjust straps by pulling adjuster until respirator fits comfortably.



4. Place fingertips at the top of nose clip. Use both hands to mold the nose clip to the shape of nose by pushing inward while moving fingertips down both sides of the nose clip.

5. Perform a User Seal Check prior to each wearing. To check the respirator-to-face seal, place both hands completely over the respirator and inhale sharply. Be careful not to disturb the position of respirator. If any leakage is detected, readjust the position of respirator according to Step #3 and #4.

Removal Instructions

Cup respirator in hand to maintain position of face and pull bottom strap over head. Still holding the respirator in position, pull top strap over head and remove respirator.

Dispose of used product in accordance with applicable regulations.