

FITTING INSTRUCTIONS:

F2
Folding Cup Roll Style Headband



STEP 1

Unfold the respirator in your hands with your fingertips on the nosepiece.



STEP 2

Holding the headbands in one hand above the respirator, grab the bottom of the respirator in the other hand.



STEP 3

Cup the bottom of the respirator under your chin, with nosepiece facing upwards.



STEP 4

Pull both headbands over your head and rest them on your neck, under your ears.



STEP 5

Pull the top headband up over your ears and rest it on the crown of your head.



STEP 6

Firmly press down the nosepiece to conform to the shape of your nose. Adjust and secure all edges and then perform a seal check.

PERFORMING A SEAL CHECK:

It is necessary to perform a seal check on the respirator every time you wear it.

Cup both hands around the respirator and inhale and exhale forcefully several times. If you notice any air leakage, adjust the nosepiece and headbands until the leakage is corrected. Once corrected, repeat the seal check again.

WARNING The respirator should not be used without a proper fit. This respirator helps protect against certain contaminants. It may not eliminate the risk of contracting disease or infection. Misuse may result in sickness or death.
Materials which may come into contact with the wearer's skin could cause allergic reactions to susceptible individuals.

Learn more at [BNX.com](https://www.bnx.com)



Manufactured by
BNX Converting, LLC.
16727 Park Row, Houston, TX 77084, USA

