

**Before using the respirator, the wearer must read and understand the User Instructions. Please keep these instructions as reference.**

**Intended Use:** This product may be used to reduce exposure to airborne infectious agents including *M. tuberculosis* and certain other airborne particles. They are designed to be fluid-resistant to protect against the splash and splatter of blood and other bodily fluids, and they are intended for use by operating room personnel and other healthcare workers.

Before occupational use of this respirator, a written respiratory protection program must be implemented meeting all the local government requirements. In the United States, employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training, and fit testing.

### User Instructions

#### Application

1. Make sure your hands are clean and dry before putting on a respirator.
2. Inspect the respirator for damage before use. If it appears damaged, dirty, or damp, do not use it.
3. Hold the respirator in one hand with the nose piece at your fingertips.
4. Place the respirator under your chin with the nose piece at the top.
5. With the other hand pull the top strap over the head, placing the strap near the crown. Then, pull the bottom strap over and place it at the back of neck below ears. Do not crisscross the straps. Make sure straps lay flat and are not twisted.
6. Place fingertips from both hands at the top of the nose piece. Press down on both sides of the nose piece to mold it to the shape of nose.

#### Seal Check

1. Place your hands on the respirator, covering as much of it as possible, then breathe out. If air leaks out from the edges of the respirator, or if glasses fog up, it is not snug. Adjust the respirator and try again.
2. Ensure your breath passes through the respirator and not around its edges.
3. If a tight seal cannot be achieved, try a different size or style.
4. This must be done every time a respirator is donned.

#### Removal

1. Without touching the respirator, lift the bottom strap from around the neck over the head. Then, lift and remove the top strap.
2. **Discard respirator after use, do not reuse respirators.**
3. Wash hands with soap and water after removing respirator.

#### Protection

N95 - Particulate Filter (95% filter efficiency level) Effective against particulate aerosols free of oil; time use restrictions may apply.

#### S - Special or Critical User Instructions

This respirator has been approved as a NIOSH N95 filtering facepiece respirator, for use in healthcare settings, as a Surgical N95 Respirator conforming to recognized standards for biocompatibility, flammability, and fluid resistance at 80 mmHg.

#### Storage Conditions and Shelf Life

Store respirators in the original packaging, away from areas with excessive moisture, dust, sunlight, chemical exposure or extreme temperatures. In original packaging at RH<80% and temperatures from -4°F(-20°C) to 86°F (30°C), the product may be used until date specified on packaging next to the hourglass icon.

