## INSTRUCTIONS / INSTRUCCIONES

**"S"- Special Instructions Section**

The RPD914P100 Air Purifying Filtering Facepiece P100 respirator has been manufactured by Venus Safety & Health Pvt. Ltd. For DSI Safety Inc. under TC-84A-7887.

**Pre-use instructions:**

- Before occupational use of this respirator, a written respiratory protection program must be implemented meeting all the local government requirements. In the United States, employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training and fit testing.
- Each time before putting on this respirator, inspect the facepiece to assure that there is no holes in the breathing zone.

**To Put On This Respirator:**

1. Hold the respirator in your hand with the metal adjusting strip away from you, allowing the headbands to hang below your hands.
2. Press the respirator firmly against your face with the narrow end on the bridge of your nose.
3. Pull the bottom strap over your head and position it below your ears. Pull the top strap over your head and position it high on the back of your head.
4. Adjust the bottom strap lengths using both hands to pull the ends evenly towards the back of your head. Straps may be loosened by pushing the back of the buckle outward.
5. Using both hands, mold the metal strip to the shape of your nose.
6. Check the facepiece-to-face seal by cupping both hands over the respirator and inhale vigorously. A negative pressure should be detected inside the facepiece. If air leaks in from around your nose, reform the metal strip. If air leaks in from around the edges, reposition the straps. Repeat as necessary until no air leak is detected, and you have a satisfactory seal. If you cannot achieve a good face-to-facepiece seal do not enter the contaminated area and see your supervisor.

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