USER'S INSTRUCTIONS

MYCO Medical Supplies, Inc.
N95 Respirator

MYTN01-11 (Medium)
MYTN01-12 (Small)

MYCO Medical Supplies, Inc.
2015 Production Drive, Apex, NC 27539
USA
TEL : +1-919-460-2535

TN01PL371UIId
[Rev. d] 06/30/2020

USER'S INSTRUCTIONS

WARNING!

Although this respirator is useful to protect the wearer from particulates, including those containing infectious disease, the risk for sickness or infection is not completely eliminated. Failure to read and follow all warnings and instruction will reduce respirator performance and may result in sickness and / or death. Before using, please read user's instruction.

[USE INSTRUCTION]

This respirator does not supply oxygen and must not be used in atmospheres containing less than 19.5% oxygen.

Use for protection against particulate aerosols free of oil.

Do not use for gasses, vapors, asbestos, paint spray or any oil-based aerosols environments.

Before occupational use of this respirator a written respiratory protection program must be implemented meeting all the local government requirements. In the United States employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training, and fit testing.

Leave the use environment immediately in the case of dizziness, nausea, breathing difficulty or other symptoms of physical distress.

Don a new respirator if a proper fit cannot be obtained, if breathing becomes difficult, or the respirator becomes damaged. Only for single use.

Store unused respirators in a clean, cool, dark, dry place.

S-[SPECIAL OR CRITICAL USER'S INSTRUCTIONS]

These MYTN01-11 and MYTN01-12 filtering face piece respirators" have been manufactured by Approval Holder San-M Package Co., Ltd. for MYCO Medical Supplies, Inc. under TC-84A-3348.

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