

# AERO PRO CO., LTD.

AP112V

## N95 Particulate Filtering Facepiece Respirator

### ■ FITTING INSTRUCTION

1.



2.



3.



4.



1. Make sure the respirator, including its head straps and nose clip, is not damaged or stained as it may negatively affect the efficiency of the respirator.
2. Place the respirator under your chin with the molded nose area up and pull the headbands over your head.
3. Position the lower head strap around the neck, below your ears, and the upper head strap around the back of your head.
4. Using both hands, adjust the nose clip to the shape of the bridge of the nose to ensure a close fit and good seal.

■ Change respirator immediately if breathing becomes difficult or respirator becomes damaged or distorted, or a proper face fit cannot be maintained. Remember, careful observance of these fitting instructions is an important step in safe respirator use.

■ To check the correct fit, cup with both hands the respirator and inhale sharply. If you feel any inward leaks around the perimeter of the respirator, stretch the headbands again by pulling them backwards. Then correct the position of the respirator and check again.

**WARNING!** If you can't achieve a proper seal, do not use the respirator.