



DLP Medical Products, Corp.
Phone: 1-(210)-858-8352
San Antonio Texas U.S.A.



1. Protection

N95 - Particulate Filter (95% filter efficiency level) effective against particulate aerosols. Free of oil. Use restrictions may apply.

This respirator is only approved in the following configuration:

| Filter | Protection | Respirator | Service and Limitation |
|--------|------------|------------|------------------------|
| 20 | NIOSH | N95C | 20-3000 |

5. Special or Critical User's Instructions: The model N95C filtering facepiece respirator has been manufactured by Approval Holder MAKRITE for Company DLP under TC-84A-3325.

J. Cautions and Limitations

- A - Not for use in atmospheres containing less than 18% oxygen.
- B - Not for use in atmospheres (immediately dangerous to life or health).
- C - Do not exceed maximum use concentration established by regulatory standards.
- J - Failure to properly use and maintain this respirator could result in injury or death.
- M - All approved respirators shall be selected, fitted, used, and maintained in accordance with NIOSH and other applicable regulations.
- N - Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer. Refer to users instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
- P - NIOSH does not evaluate respirators for use in surgical masks.
- S - Special or critical user's instructions and specific use limitations apply. Refer to the instructions before donning.

Before use, must be fit tested either quantitatively or qualitatively by a qualified fit tester before using the respirator in a contaminated environment. Refer to the OSHA guidelines, 29 CFR.

Service Time Limitation

1. Do not use the same respirator for more than 8 hours, continuous or intermittent, in all workplaces that could result in high filter loading.
2. Leave contaminated areas and remove respirator. If breathing becomes difficult proper fit cannot be obtained, or if respirator is damaged, soiled, or distorted use the same respirator again once removed.
3. The maximum service time of this respirator is five (5) consecutive calendar days (including days of non-use) beginning from the first day of use.

DLP Medical Products, Corp.
Phone: 1-(210)-858-8352
San Antonio Texas U.S.A.



CAUTION - WARNINGS - PRECAUTIONS - ADVISORY INFORMATION

This limited warranty is made in lieu of the expression of workmanship, fitness for particular purposes and all other warranties, express or implied. There are no other warranties which extend beyond the description on the face hereof. The physical structure and configuration of it will be not be protected.

Exclusive Remedies: Damages for the breach of this limited warranty are limited to the replacement of such quantity of the products proved to be defective. Damages, including incidental, special, or consequential, arising out of such, any, or delay, or the inability to use products by the user.

HOW TO USE:

1. Read and read fitting instruction and use instructions before use.
2. This respirator does not protect against the risk of contracting disease or infection.
3. Before use, the user must be trained in the proper use to accordance with applicable safety and health standards for the contained and exposed level in the work area.



Fitting Instructions:

- step 1**
Hold respirator in your hand as shown, with the nosepiece at your fingertips and the headbands hanging freely below your hand.
- step 2**
Place respirator under your chin with the nosepiece up. Pull the top headband to top back of head.
- step 3**
Pull the bottom headband overhead and place around neck below ears.
- step 4**
Adjust headbands for comfortable fit. Push the nosepiece inward with your fingertips to conform to the shape of your nose.



Caution:
The respirator should be fit checked for each wearing. Place both hands completely over the respirator and exhale. If air leaks around nose or the edges, adjust the nosepiece and/or headbands until a good fit is achieved.

Filtering facepieces are to be inspected prior to each use assure there are no holes in the filtering media other than the pores around staples and no damage has occurred. Exchanged body resulting from ripped or torn filter material around staples are considered damage.



DLP Medical Products, Corp. Made in China

LOT NO.