

Fido disposable particulate respirator



N95



FILTERING FACEPIECE
APPROVAL LABEL



FIDO MASK CO., LTD.
No. 8-308-4th Fl. Dong 41 Road,
Miaoli City 305, Taiwan
TEL: 886-3-571-1171

FIDO RESPIRATOR APPROVAL LABEL IN THE FOLLOWING CONFIGURATION

FC	TradeMark	Respirator	Caution and Limitation ¹⁾
88A-22750	N95	FIDO	DISPOSABLE

1. Precaution

NIOSH (National Institute for Occupational Safety and Health) is the official body for filtering respirator approval. Please refer to the instructions for use for more detailed information.

2. Caution and Limitation

1. Not for use in environments containing less than 10% oxygen.
2. Not for use in atmospheres immediately dangerous to life or health.
3. Do not exceed maximum use concentration established by regulatory standards.
4. Subject to proper use and inspection, the product will last for a single work shift.
5. All approved respirators shall be selected, fitted, used, stored, and discarded in accordance with OSHA, OSHA's and other applicable regulations.
6. Not for use in atmospheres, aerosols, or vapors. Use only over respiratory protection that is designed to specifically filter aerosols.
7. Refer to user's instruction, and the manufacturer's information on the front and back of the respirator, for use and storage instructions.
8. NIOSH does not evaluate respirators for use in unapproved media.

F750 N95 Particulate Respirator

User's Instructions

WARNING TO USER

Observe the safety instructions, warnings and labels on the F750 N95 Particulate Respirator and the instructions. This respirator and the filter are not intended to be used as a substitute for a complete safety program or to provide protection from environmental conditions that are not covered by the approved safety program. The respirator does not protect against chemical, biological, radioactive or ionizing radiation, or protection from conditions that require special eye protection. Failure to follow the instructions, warnings and labels, combined with the use of this respirator, may result in serious injury or death.

The F750 N95 Particulate Respirator is a Class II, Type I, N95 particulate respirator with a lead and gasket face seal and an exhalation valve.

DO NOT USE FOR:

The F750 N95 respirator should not be used if you require additional information on the use of environmental control (ECM) respirators, such as high efficiency particulate respirators, or for other uses not specifically approved by NIOSH.

USE PRECAUTIONS:

1. Only use respirators for the purpose for which they were designed. Do not use for purposes other than those for which they were designed. Do not use for purposes for which they were not designed. Do not use for purposes for which they were not designed. Do not use for purposes for which they were not designed.
2. Do not use the respirator in atmospheres containing less than 10% oxygen.
3. Do not use the respirator in atmospheres immediately dangerous to life or health.
4. Do not use the respirator in atmospheres containing less than 10% oxygen.
5. Do not use the respirator in atmospheres immediately dangerous to life or health.
6. Do not use the respirator in atmospheres containing less than 10% oxygen.

USE LIMITATIONS:

1. This respirator has not been approved for use in atmospheres containing less than 10% oxygen.
2. Do not use the respirator in atmospheres immediately dangerous to life or health.
3. Do not use the respirator in atmospheres containing less than 10% oxygen.
4. Do not use the respirator in atmospheres immediately dangerous to life or health.
5. Do not use the respirator in atmospheres containing less than 10% oxygen.
6. Do not use the respirator in atmospheres immediately dangerous to life or health.

STORAGE:

1. Do not store respirators in their original packaging in a clean, dry, cool environment that is not subject to extreme temperatures, excessive humidity, dust, or chemicals.
2. Do not store respirators in their original packaging in a clean, dry, cool environment that is not subject to extreme temperatures, excessive humidity, dust, or chemicals.
3. Do not store respirators in their original packaging in a clean, dry, cool environment that is not subject to extreme temperatures, excessive humidity, dust, or chemicals.

If your respirator is damaged, it should not be used for the purpose for which it was designed.



Fido disposable particulate respirator



Fitting Instructions:

Before using this respirator, carefully inspect the respirator to ensure it meets the requirements for use in an approved and functioning condition. Remove the entire respirator if you find any of the following problems:

1. If the fit of the respirator to your face is not secure, do not use it until you have corrected the fit.
2. If the respirator does not fit snugly against the face, do not use it until you have corrected the fit.
3. If the respirator does not fit snugly against the face, do not use it until you have corrected the fit.
4. If the respirator does not fit snugly against the face, do not use it until you have corrected the fit.

Note: Carefully inspect the respirator before use. Do not use if you find any of the following problems:

1. If the fit of the respirator to your face is not secure, do not use it until you have corrected the fit.
2. If the respirator does not fit snugly against the face, do not use it until you have corrected the fit.
3. If the respirator does not fit snugly against the face, do not use it until you have corrected the fit.
4. If the respirator does not fit snugly against the face, do not use it until you have corrected the fit.

For more information, visit the NIOSH website at www.niosh.gov.

CAUTION: Do not use.