

Willson®

A Bacou-Dalloz Company
3001 South Susan Street
Santa Ana, CA 92704

1-800-821-7236

DESCRIPTION

The ONE-Fit HC-NB095 particle respirator and surgical mask is a class N95 disposable mask designed to provide filter efficiency of 95% or greater against particulate aerosols free of oil. The respirator provides greater than 99% Bacterial Filtration Efficiency (BFE) against micro-organisms generated by a wearer, and meets CDC guidelines for M. tuberculosis exposure control. It can fit a wide range of face sizes, is fluid resistant and suitable for surgical environments.

Note: Bacterial Filtration Efficiency (BFE) tested in accordance with MIL M36954C, ASTM F2101. Fluid resistance conforms to ASTM F1862, 160 mmHg.

INTENDED USE

The Willson ONE-Fit HC-NB095 Healthcare Particulate Respirator and Surgical Mask is a NIOSH-approved N95 single use respirator intended for use by health care personnel during medical/surgical procedures to protect both the wearer and the patient by protecting the wearer against the spatter of blood and other potentially infectious materials and reducing the transfer of micro-organisms and other airborne particulate matter.

The Willson ONE-Fit also meets the CDC guidelines for TB exposure control within healthcare facilities and is intended for use as an isolation mask.



FITTING INSTRUCTIONS

Follow these instructions each time respirator is worn.

1. Cup the respirator in the hand, allowing the head straps to hang freely.
2. Hold the respirator under the chin with the nosepiece facing outwards.
3. Place the lower head strap around the neck below the ears.
- 4/5. Holding the respirator against the face with one hand, place the top head strap above the ears, around the crown of the head.
6. Place hands on each side of the respirator and move slightly right, left, up and down, to adjust the position of the respirator and achieve the most optimal fit on the face.

FIT-CHECK THE FACE-SEAL AS FOLLOWS:

1. Place both hands over the respirator - do not disturb the position of the respirator;
2. Inhale sharply (a negative pressure should be felt inside the respirator);
3. If you detect air-leaks, readjust the head straps and/or the nose-piece.

WARNING: If you cannot achieve satisfactory fit, do NOT enter the contaminated area. See your Supervisor.

USE INSTRUCTIONS:

1. Failure to follow all instructions and limitations concerning the use of the respirator and/or failure to wear the respirator for the whole duration of exposure to contaminants can seriously reduce the performance characteristics of the respirator and lead to illness, injury or death.
2. Before use, the wearer must first be trained by the employer in the correct use of the respirator in accordance with applicable safety and health standards.
3. OSHA standards including 29CFR1910.134(e)(5), (f) and CSA standard Z94.4-93 require that the wearer be fit-tested.
4. Leave the contaminated area if breathing becomes difficult or if the respirator becomes damaged.
5. Discard the respirator and replace with a new one if:
 - excessive clogging of the respirator causes breathing difficulty.
 - the respirator becomes damaged.
 - the respirator becomes contaminated with blood or bodily fluids
6. Discard the respirator after every use when used for surgical procedures.
7. Do NOT use the respirator with beards or other facial hair that interferes with direct contact between the face and the edge of the respirator, or any other conditions that may prevent a good face-seal.
8. This respirator is NOT to be used on children.

WARNING:

This respirator helps protect against certain contaminants. It may not eliminate the risk of contracting disease or infection. **Misuse may result in sickness or death.**

The straps of this respirator are made from latex-free rubber. However, individuals who are highly sensitive to natural rubber latex may have an allergic reaction.

FOR MORE INFORMATION or assistance on Willson® products, contact your local Willson distributor or Bacou-Dalloz.