

四色



130 MM

N95

LATEX
FREE



20 Respirators



ARAGON PTY LTD



National Institute for
Occupational Safety and Health
NIOSH

This respirator is approved only in the following configuration:

TC =	Protection ¹	Respirator	Cautions and Limitations ²
86A-5111		ARAGON N95	
	N95	Y	ABC INADRS

Protection

N95 - Particulate Filter (95% filter efficiency level) effective against particulate aerosols free of oil; lime use restrictions may apply.

2. Cautions and Limitations

- A - Not for use in atmospheres containing less than 19.5% oxygen.
- B - Not for use in atmospheres immediately dangerous to life or health.
- C - Do not exceed maximum use concentrations established by regulatory standards.
- J - Failure to properly use and maintain this product could result in injury or death.
- M - All approved respiratory use: to be isolated, fitted, used, and maintained in accordance with MSDA, OSHA and other applicable regulatory requirements.
- N - Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- O - Refer to user's instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
- P - NIOSH does not evaluate respirators for use as surgical masks.
- S - Special or critical User's Instruction and/or specific use limitations apply.

Service time limitations

1. Do not use the same respirator for more than 8 hours, continuous or intermittent, in dirty workplaces that could result in high fiber loading.
2. Leave contaminated area and remove respirator if breathing becomes difficult or proper fit cannot be obtained or the respirator is damaged, soiled, or distorted. Do not use the same respirator again once removed.
3. The maximum service time of this respirator is five (5) consecutive calendar days (including days of non-use), beginning from the first day of use.



LIMITED WARRANTY

LIMITED WARRANTY
IMPORTANT NOTICE TO PURCHASER

This limited warranty is made in lieu of the warranties of merchantability, fitness for particular purposes and all other warranties, express or implied. There are no other warranties which extend beyond the description on the face hereof. The physical standards and specifications of Makrite will be met by products sold.

Exclusive Remedies: damages for the breach of this limited warranty are limited to the replacement of such quantity of Makrite, products proved to be defectively manufactured. Except as provided above, Makrite shall not be liable or responsible for any loss, damage, or liability, direct, indirect, incidental, special, or consequential, arising out of sale, use, or misuse, or the inability to use products by user.

ARAGON PTY LTD
www.aragon-mine.com.au
www.aragon.com.au
PO Box 1009, (Doodlakine South, NT) 2744-Australia

Made in China

WARNINGS

2. The respirator does not protect against the risk of contracting disease or infection.
3. Before occupational use of this respirator, a written respiratory protection program must be implemented meeting all the local government requirements. In the United States employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training and fit testing.



STEP 1

STEP 1 Hold the respirator in hand with the nosepiece at your fingertip, allowing the headbands to hang freely below your hand.



STEP 3

Press the respirator firmly against your face with the nosepiece on the bridge of your nose.



STEP 3

Stretch and position the top band high on the back of the head. Stretch the bottom band over the head and position below your ears.



000000

STEP 4
Using both hands, mold nosepiece to the shape of your nose.



STEP 5

To test fit: a) Cup both hands over the respirator being careful not to disturb position, and b) inhale vigorously. If air leaks around the edges, reposition the straps or adjust strap tension for better fit.

Caution: The respirator should be fit checked prior to each wearing.

ARAGOON PTY LTD
 www.aragoon.com.au
 info@aragoon.com.au
 6/1, Bldg. 4/100, Southport Green, 5008, 21/05/2010

LOT #

Expiration Date:

