## INSTRUCCIONES / INSTRUCTIONS / INSTRUCTIONS

**Español**

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**English**

Before the occupational use of this respirator, a written respiratory protection program must be implemented that meets all the requirements of Regulation 522-06 of the Government of the Dominican Republic. In the United States, employees must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training and fit testing. The user of this respirator must be properly trained by the supervisor as to the correct use of the respirator and according to the corresponding safety and hygiene standards. The respirator should not be used for more than 8 hours.

Fitting Instructions:

1. Hold the respirator cup in your hand, allowing the elastics to hang freely.
2. Place the respirator under the chin with the nose piece facing outwards.
3. Place the lower head strap around the neck below the ears.
4. Holding the respirator against the face with one hand, place the top head strap above the ears, around the crown of the head.
5. Place your hands on each side of the respirator and gently move it right, left, and up, down to achieve optimum respirator fit. Fit Test: Exhale sharply, this will produce positive pressure inside the respirator causing it to expand. If air flows around the nose, tighten the metal strip. If air flows around the edges, replace the bands for a better fit. Repeat as necessary until no air leak is detected and has a satisfactory adjustment. If you CAN NOT achieve a proper fit, DO NOT enter the contaminated area. Consult your supervisor for proper respiratory protection.

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