

# RESPIRATOR FITTING INSTRUCTION:

*Please carefully follow these fitting instructions during each use to achieve proper fit.*



**1.** Hold the respirator in hand with the nosepiece at your fingertips, allowing the elastic straps to hang freely below your hand.



**2.** Press the respirator firmly against your face with the nosepiece on the bridge of your nose.



**3.** Stretch and position the top band high on the back of the head. Stretch the bottom band over the head and position below your ears.



**4.** Using both hands, mold nosepiece to the shape of your nose.



**5.** To test fit: a) Cup both hands over the respirator being careful not to disturb position, and b) inhale vigorously. If air leaks around the edges, reposition the straps.