



Particulate Respirator

QB-920-N95/QB-920S-N95

Particulate Respirator
Recommended for grinding, sanding, sawing, insulating and sweeping.

QB-920-N95/QB-920S-N95

WARNING

This respirator can help reduce inhalation exposures to certain airborne biological particles but cannot eliminate the risk of contracting infection, illness or disease.**Misuse may result in sickness or death.** For proper use, see supervisor, or *User Instructions*, or write to :

QUEBEE SAFETY TECHNOLOGIES (INT’L) PTE LTD
239 Ubi Ave 4
Intrepid Warehouse Complex
Singapore 408821
Tel : (65) 6347 1689

IMPORTANT

Before use, the wearer must read and understand these User Instructions.

Description

The QB-920-N95 / QB-920S-N95 Health Care Particulate Respirator is designed to help provide respiratory protection for the wearer. This product has a filter efficiency level of 95% or greater against particulate aerosols free of oil. It is fluid resistant, disposable and may be worn in surgery. It can fit a wide range of face sizes.

1 Tested against a 0.3 micron particle (count median diameter) per 42 CFR 84.

INTENDED USE

Particles such as those from grinding, sanding, sawing, bagging, insulating, sweeping, or processing minerals, iron ore, coal, flour, wood, metal, pollen, and certain other substances. Liquid or non-oil based particles from sprays that do not also emit oil aerosols or vapors.

NOT FOR USE

Gases and vapors, oil aerosols, asbestos, arsenic, cadmium, lead, 4,4'-methylene dianiline (MDA) or sandblasting.

TIME USE LIMITATION

If respirator becomes soiled, damaged, or breathing becomes difficult, leave the contaminated area immediately and replace the respirator.

USE INSTRUCTIONS

- Before occupational use of this respirator a written respiratory protection program must be implemented meeting all the local government requirements. In the United States employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training, and fit testing.
- Failure to follow all instructions and limitations on the use of this respirator and/or failure to wear this respirator during all times of exposure can reduce respirator effectiveness, and may result in sickness or death.
- Filtering facepieces are to be inspected prior to each use to assure there are no holes in the breathing zone other than the punctures around staples and no damage has occurred. Enlarged holes resulting from ripped or torn filter material around staple punctures are considered damage. Immediately replace respirator if damaged. Staple perforations do not affect NIOSH approval.
- Store the respirator away from contaminated areas when not in use.
- Leave the contaminated area immediately and contact supervisor if dizziness, irritation, or other distress occurs.
- The particles which can be dangerous to your health include those so small that you cannot see them.
- Dispose of used product in accordance with applicable regulations.

LIMITATIONS

- This respirator does not supply oxygen. Do not use in atmospheres containing less than 19.5% oxygen.
- This respirator is designed for occupational/professional use by adults who are properly trained in their use and limitations. Respirators are not designed to be used by children.
- Do not alter, abuse or misuse this respirator.
- Do not use when concentrations of contaminants are immediately dangerous to life and health, are unknown or when concentrations exceed 10 times the permissible exposure limit (PEL) or according to specific OSHA standards or applicable government regulations, whichever is lower.
- Do not use with beards or other facial hair or other conditions that prevent a good seal between the face and the sealing surface of the respirator.
- Respirators can help protect your lungs against certain airborne contaminants. However, they will not prevent entry through other routes such as the skin, which would require additional personal protective equipment (PPE).
- Individuals with a compromised respiratory system, such as asthma or emphysema, should consult a physician and complete a medical evaluation prior to use.



FITTING INSTRUCTIONS

Must be followed each time respirator is worn.

- Cup the respirator in your hand, with the nosepiece at your fingertips, allowing the headbands to hang freely below your hand.
- Position the respirator under your chin with your nosepiece up. Pull the top strap over your head resting it high at the top back of your head. Pull the bottom strap over your head and position it around the neck below the ears.
- Place your fingertips from both hands at the top of the metal nosepiece. Using two hands mold the nose area to the shape of your nose by pushing inward while moving your fingertips down both sides of the nosepiece.

▲ Pinching the nosepiece using one hand may result in improper fit and less effective respirator performance. Use two hands.

REMOVAL INSTRUCTIONS

- Cup respirator in hand to maintain position on face.
- Pull bottom strap over head.
- Still holding respirator in position, pull top strap over head and remove respirator.

NIOSH Approved: N95

At least 95% filtration efficiency against solid and liquid aerosols that do not contain oil.

QUEBEE SAFETY TECHNOLOGIES (INT’L) PTE LTD
239 Ubi Ave 4
Intrepid Warehouse Complex
Singapore 408821
Tel : (65) 6347-1689

NIOSH
National Institute for Occupational Safety and Health

THIS RESPIRATOR IS APPROVED ONLY IN THE FOLLOWING CONFIGURATION:

TC-	Protection ¹	Respirator	Cautions and Limitations ²
		QB-920-N95	
84A-7867	N95	X	ABCJMNOPS

1. PROTECTION

N95 - Particulate Filter (95% filter efficiency level) effective against particulate aerosols free of oil ; time use restrictions may apply
--

2. CAUTIONS AND LIMITATIONS

- A - Not for use in atmospheres containing less than 19.5 percent oxygen.
- B - Not for use in atmospheres immediately dangerous to life or health.
- C - Do not exceed maximum use concentrations established by regulatory standards.
- J - Failure to properly use and maintain this product could result in injury or death.
- M - All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA,OSHA, and other applicable regulations.
- N - Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- O - Refer to user's instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
- P - NIOSH does not evaluate respirators for use as surgical masks.
- S - Special or critical User's Instructions and/or specific use limitation apply. Refer to User's Instructions before donning.

QUEBEE SAFETY TECHNOLOGIES (INT’L) PTE LTD
239 Ubi Ave 4
Intrepid Warehouse Complex
Singapore 408821
Tel : (65) 6347-1689

NIOSH
National Institute for Occupational Safety and Health

THIS RESPIRATOR IS APPROVED ONLY IN THE FOLLOWING CONFIGURATION:

TC-	Protection ¹	Respirator	Cautions and Limitations ²
		QB-920S-N95	
84A-7871	N95	X	ABCJMNOPS

1. PROTECTION

N95 - Particulate Filter (95% filter efficiency level) effective against particulate aerosols free of oil ; time use restrictions may apply
--

2. CAUTIONS AND LIMITATIONS

- A - Not for use in atmospheres containing less than 19.5 percent oxygen.
- B - Not for use in atmospheres immediately dangerous to life or health.
- C - Do not exceed maximum use concentrations established by regulatory standards.
- J - Failure to properly use and maintain this product could result in injury or death.
- M - All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA,OSHA, and other applicable regulations.
- N - Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- O - Refer to user's instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
- P - NIOSH does not evaluate respirators for use as surgical masks.
- S - Special or critical User's Instructions and/or specific use limitation apply. Refer to User's Instructions before donning.

S - Special or Critical User's Instructions: The model QB-920-N95 & QB-920S-N95 filtering facepiece respirator has been manufactured by Approval Holder Makrite for Company QUEBEE SAFETY TECHNOLOGIES (INT’L) PTE LTD under TC-84A-7867 & TC-84A-7871.