

Medicom SafeMask Architect Pro Respirator User Instructions

Models 203214, 203314 and 203414

Medicom Surgical N95 Respirators for single use meet the requirements of the 42 CFR Part 84 standard. These filtering facepiece respirators are recommended for reduced exposure against particulate aerosols free from oil. They have also been evaluated for use in healthcare settings according to recognized standards for biocompatibility (ISO 10993), flammability (16 CFR 1610), and fluid resistance (as per ASTM F1862 120 mm Hg).

DONNING & FITTING INSTRUCTIONS

Follow these instructions each time respirator is worn.



Pull open the top and bottom edges of the respirator to fully open.



Bend the nose clip into a slight "V curve" so that the respirator forms a cup.



Hold the single split comfort band and pull apart to create two separate comfort bands and rest these atop the respirator.



Position the bottom edge of the respirator under the chin and the top edge on the nose. Bring the two comfort bands over the head.



Position one comfort band across the crown of the head.

Medicom SafeMask Architect Pro Respirator User Instructions

Models 203214, 203314 and 203414



Position the other comfort band around the back of the neck.



Mold the nose clip to the shape of the face, running the fingertips of both hands from the top of the nose clip down both sides while pressing inward. Always use both hands.



FIT-CHECK THE FACE-SEAL AS FOLLOWS: Place both hands over the respirator - do not disturb the position of the respirator. Exhale sharply (a positive pressure should be felt inside the respirator). If you detect air-leaks around edges of the respirator, readjust the elastic bands and/or the nose clip.

FIT-CHECK THE FACE-SEAL AS FOLLOWS:

Place both hands over the respirator - do not disturb the position of the respirator. Exhale sharply (a positive pressure should be felt inside the respirator). If you detect air-leaks around edges of the respirator, readjust the elastic bands and/or the nose clip.

WARNING: If you cannot achieve satisfactory fit, do NOT enter the contaminated area. See your Supervisor. This respirator reduces exposure against certain contaminants. It may not eliminate the risk of contracting disease or infection. Misuse may result in sickness or death.

DOFFING INSTRUCTIONS

Remove both comfort bands by bringing them back over the head and remove the respirator from the face.

STORAGE INSTRUCTIONS

When not in use store the respirator in its original packaging or a plastic bag in a clean and dry environment.

AMD Medicom, Inc.
Montréal, Québec, Canada
1-800-361-2862

AP203-1UI Rev 00 October 2021

Medicom SafeMask Architect Pro Respirator User Instructions

Models 203214, 203314 and 203414

USER INSTRUCTIONS AND USE LIMITATIONS

1. Failure to follow all instructions and limitations concerning the use of the respirator and/or failure to wear the respirator for the whole duration of exposure to contaminants can seriously reduce the performance characteristics of the respirator and lead to illness or injury or death.
2. Before use, the wearer must first be trained by the employer in the correct use of the respirator in accordance with applicable safety and health standards.
3. Before use, ensure that the respirator facepiece and the elastic bands are not damaged.
4. The respirator does not require any maintenance. It has been designed for one single continuous work session, i.e. for 8 hours as a maximum.
5. Before occupational use of this respirator, a written respiratory protection program must be implemented meeting all the local government requirements. In the United States, employers must comply with OSHA 29CFR1910.134, which include medical evaluation, training, and fit testing. In Canada, employers must comply with CSA standard Z94.4-93.
6. Discard the respirator and replace with a new one if:
 - a. excessive clogging of the respirator causes breathing difficulty.
 - b. the respirator becomes damaged.
7. Leave the contaminated area if breathing becomes difficult; dizziness, irritation or other distress occurs; you smell, taste or sense irritation from the contaminants; or the respirator becomes damaged.

Do NOT use the respirator or enter or stay in a contaminated area under the following circumstances:

- atmosphere contains oil aerosols
- for protection against gases or vapors

Do NOT modify or misuse the respirator. Do NOT use the respirator with beards or other facial hair that interferes with direct contact between the face and the edge of the respirator, or any other conditions that may prevent a good face-seal.

AMD Medicom, Inc.
Montréal, Québec, Canada
1-800-361-2862

AP203-1UI Rev 00 October 2021