

Fido disposable particulate respirator



N95



FILTERING FACEPIECE  
APPROVAL LABEL



FIDO MASKS CO., LTD.  
No. 8-5 She-Ling Li, Yuan-Li Town,  
Miao-Li County 35872, Taiwan  
TEL: 886-37-741471

THIS RESPIRATOR IS APPROVED ONLY IN THE FOLLOWING CONFIGURATION:

TC-	Protection <sup>1</sup>	Respirator	Cautions and Limitations <sup>2</sup>
84A-FITF750	N95	F750 X	ABCDMNOP

#### 1. Protection

N95 - Particulate Filter (95% filter efficiency level)  
Effective against particulate aerosols free of oil.  
Time use restrictions may apply.

#### 2. Cautions and Limitations

- A - Not for use in atmospheres containing less than 19.5% oxygen.
- B - Not for use in atmospheres immediately dangerous to life or health.
- C - Do not exceed maximum use concentrations established by regulatory standards.
- J - Failure to properly use and maintain this product could result in injury or death.
- M - All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA and other applicable regulations.
- N - Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- O - Refer to User's Instructions, and/or maintenance materials for information on use and maintenance of these respirators.
- P - NIOSH does not evaluate respirators for use as surgical masks.

## F750 N95 Particulate Respirator

### User's Instructions

#### WARNING TO USER

Before use, the user must read, understand and follow all User's Instructions, warnings and the limitations of this respirator. This respirator provides limited protection and must be properly selected and fit tested by a competent person to provide protection from known contaminants according to OSHA 29 CFR 1910.134 or applicable safety and health standards. The respirator's face seal must be checked each time before entering a contaminated area and the respirator must be worn during all times of exposure to contaminants. Failure to follow all of these instructions, warnings and limitations, contained on this package may reduce the effectiveness of this respirator and result in serious injuries or DEATH!

The F750 N95 respirator with a filtration efficiency of at least 95% against solid or liquid particles that do not contain oil.

#### DO NOT USE FOR:

The F750 N95 respirator should not be used for gases, vapors, oil aerosols or when a concentration of contaminants exceeds OSHA or applicable safety and health standards' maximum use concentrations, are unknown or immediately dangerous to life and health.

#### USE INSTRUCTIONS:

- Before occupational use of this respirator a written respiratory protection program must be implemented meeting all the local government's requirements. In the United States employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training, and fit testing.
- Prior to each use, carefully inspect the entire respirator for correct assembly and for tears or damage including the straps, filter and exhalation valve (if equipped with one). **Discard the entire respirator if any part does not function properly and get a new respirator.**
- If the respirator becomes damaged, distorted, soiled, or breathing becomes difficult or a proper respiratory face seal cannot be maintained, **leave the contaminated area immediately and get a new respirator.**
- If distress, irritation, nausea or other distress occurs, **leave the contaminated area immediately and** contact your supervisor or medical professional.
- Do not talk, smoke or eat while using this respirator.
- Properly dispose of used respirators in accordance with applicable regulations.

#### USE LIMITATIONS:

- This respirator does not supply oxygen and must not be used in atmospheres containing less than 19.5% oxygen.
- This respirator is not to be used when the user has breath, facial hair, sores or other conditions that prevent a good seal between the face and the respirator.
- This respirator is designed for occupational or professional use by adults who are properly trained in their use and limitations. **Respirators must not be used by children.**
- Particulates and other contaminants which can be dangerous to your health include those that you cannot see, taste or smell. Respirators can help to protect your lungs against certain airborne contaminants, but cannot prevent entry through other routes such as the skin. Additional types of personal protective equipment may be required.

#### STORAGE:

- Store unused respirators in their sealed, original packaging in a clean, dry, non-contaminated area that is protected from extreme temperatures, moisture, sunlight, dust and chemicals.
- Never store in a clean and reuse this respirator.
- Used respirators should not be stored and should be disposed of properly.



See information supplied by the manufacturer.



End of Shelf Life: 3 years.



Temperature range of storage conditions: -10°C to +40°C.



Maximum relative humidity of storage conditions: <70%.

Fido disposable particulate respirator



N95



#### FITTING INSTRUCTIONS:

Before fitting the respirator, carefully inspect the entire respirator to ensure that the respirator is not damaged and is functioning correctly. **Discard the entire respirator if any part does not function properly.**

- Hold the respirator in your hand with the nosepiece at your fingertips and allowing the straps to fall below the hand.
- Place respirator under the chin with the nosepiece (narrow end) up. Pull the continuous bottom strap over your head and position it around the neck below the ears.
- Raise the continuous top strap over your head and place it high to the back of the head.
- Using both hands, push the nosepiece towards your fingertips to mold and shape the nose piece to the face.
- A seal check MUST be performed each time the user enters a contaminated area. Lightly cup both hands over the front of the respirator and INHALE SHARPLY. A negative pressure should be felt inside of the respirator. If any air flows is detected around the edges of the respirator, readjust the straps. Repeat seal check until a proper seal is achieved; otherwise, see your supervisor. **NEVER enter a contaminated area without a proper seal.**

Note: Plucking the nosepiece using one hand may result in improper fit and less effective respirator performance. Use two hands.

#### Manufacturer Address:

Fido Mask Co., Ltd.  
8-5 She-Ling Li, Yuan-Li Town, Miao-Li County 35872, Taiwan

Production date: YYYY/MM

Lot No. XXXXXX

(The Lot No. XXXXXX will be as same as abbreviated Label indicated, the "XXXXXX" should be printed as purchase order number.)

Contents: 20 pcs.