

User's Instructions

6/APR/2010

Thank you for purchasing our product.

Before use, wearer must read and understand these User's instructions. Keep these User's instructions so you can see at anytime. For More information, feel free to contact us

.(☎ +82-31-429-2042~3)

Filtering half mask to protect against particles, X-100



FiNEtech Co., Ltd

B103 Hyundai I Valley, #14-1, Dang-dong, Gunpo-city, Gyeonggi-do, Korea
Tel : +82-31-429-2042~3 Fax : +82-31-429-2044

1. Usage

This product is a P1 grade filtering half mask to protect against particles.



DANGER - Misuse result in sickness or death. Do not use this product beyond usage. People who use this respirator first must be instructed by supervisor knowing how to use this. Never use in the following cases.

- in atmospheres containing less than 19.5 percent oxygen
- in case of the contaminants are unclear
- in atmospheres existing poisonous gas
- in atmospheres can be exposed to radioactive material or dioxin
- in atmospheres the contaminants are immediately dangerous to life or health



WARNING – Misuse may result in sickness or death.

- Check this respirator before use.
- Do not use this respirator found to have defects before use.
- Do not disassemble or alter this respirator.
- Do not use non-genuine parts.
- Do not use with towels or other conditions that prevent a good seal between the face and the sealing surface of the respirator.
- Replace the respirator in a safe place if breathing becomes difficult.



CAUTION – Misuse may result in property damage or injury.

- Do not take off the respirator during work.
- This respirator helps protect against certain particles.
- Do not reuse after cleaning.
- Use the respirator after checking the adhesion to the face

2. Characteristics

- 5 layer filter media for high efficiency
- Everyone can wear this respirator regardless of the face size because head string can be fixed to the face.
- In spite of a long use, this respirator is wearable.

3. Fitting instructions



- ① Adhere to the face closely by pulling the right hook and left head strings.
- ② Hang the left head strings to the right hook.
- ③ Push the nose clip with fingers of both hands to prevent the inflow of outside air.
- ④ Adjust the head strings not to feel pressure on the face.
- ⑤ Check again whether the outside air flows into the respirator when inhaling.

4. Maintenance and repair

- 1) Replacement – This respirator can't be replaced because it is disposable.
- 2) Cleaning - This respirator can't be cleaned because it is disposable.
- 3) Storage



Keep this respirator in a dry condition in order not to be transformed.

- 4) Repair - This respirator can't be repaired because it is disposable.
- 5) Discard and Recycle - Dispose of used respirator as a industrial waste. This respirator can't be recycled

Date of manufacture : Displayed separately

Seller :Finetech Co., Ltd.

Consumer counseling :+82-31-429-2042~3



FILTERING FACEPIECE
APPROVAL LABEL

FINetech Co., Ltd
B103, Hyundai I Valley, #14-1, Dang-dong,
Gunpo-city, Gyeonggi-do, 435-010, Republic of Korea
TEL : 82-(0)31- 429 – 2042 FAX : 82-(0)31- 429 - 2044

THIS RESPIRATOR IS APPROVED ONLY IN THE FOLLOWING CONFIGURATION:

TC-	Protection ¹	Respirator	Cautions and Limitations ²
		X-100	
84A-5461	N95	X	ABCJMNOP

1. Protection

N95 - Particulate Filter (95% filter efficiency level)
Effective against particulate aerosols free of oil;
time use restrictions may apply

2. Cautions and Limitations

- A - Not for use in atmospheres containing less than 19.5% oxygen.
- B - Not for use in atmospheres immediately dangerous to life or health.
- C - Do not exceed maximum use concentrations established by regulatory standards.
- J - Failure to properly use and maintain this product could result in injury or death.
- M - All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA and other applicable regulations.
- N - Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- O - Refer to users instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
- P - NIOSH does not evaluate respirators for use as surgical masks.