Respirator Use Instructions:

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84A-7853

1. Hold the respirator in hand with nosepiece at your fingertips, allowing the headbands to hang freely below your hand.
2. Press the respirator firmly against your face with the nosepiece on the bridge of your nose.
3. Stretch and position the top band high on the back of your head. Stretch the bottom band over the head and position it below your ears.
4. Using both hands, mold the metal nosepiece to the shape of your nose.
5. To test fit:
 - a) Cup both hands over the front of the respirator, being careful not to disturb its position.
 - b) Inhale sharply. If leakage is suspected, adjust position of respirator and/or strap and repeat procedure.
6. Remember, careful observance of these fitting instructions is an important step in safe respirator use.

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NOTE: Before occupational use of this respirator, a written respiratory protection program must be implemented meeting all the local government requirements. In the United States, employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training, and fit testing.

* Gerson recommended for relief against the irritating effects of nuisance levels of organic vapor or acid gas. This product does not provide any protection from organic vapors or acid gases. Nuisance level refers to concentrations not exceeding OSHA PEL or applicable government occupational exposure limits, whichever is lower.

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