

JSP 821, 822, 823

Donning Procedures

1. Hold the respirator in one hand with nosepiece at your fingertips, allowing the headbands to hang below your hand. The harness clip should be unsecured.
2. Press the respirator against your face with the nose part of the mask resting on the bridge of your nose.
3. Place the top band high at the back of your head.
4. Pull the two parts of the lower band back behind the neck and secure the clip. Adjust the tension of the strap by pulling the top band up from the head and sliding the tension adjuster down to the head to achieve a comfortable fit.
5. Test the fit of the respirator. Cup both hands over the respirator and exhale vigorously. If air leaks around the edges of the mask, tighten and/or reposition the bands for a better fit. Repeat this test.
6. Remember, careful observance of these fitting instructions is an important step for safe respirator use. If you CANNOT achieve a proper fit, DO NOT enter the contaminated area. See your supervisor for the instruction to obtain the proper fitting of your respirator.

