KOSTO

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RMK010

Approved
Homologué

Attention:
Before use, read the recommendations,
limitations as well as the application
instructions indicated on the box.

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If you CANNOT achieve proper fit, DO NOT enter the contaminated area.
See a qualified supervisor for proper fit of respiratory protection.

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KOSTO

N95 NIOSH Approved
RMK010
10 Disposable Respirators

Produits Kosto / Kosto Products
60, Grande-Duvalville, Brossard, Québec, Canada J3C 6A3
1866-720-9292

NIOSH
N95
K
DECLASSIFIED

This respirator is only approved for the following configurations:

TC	Filter/Cartridge	Model	Color and Markings
844-3004	N95	RMK010	X

1. Protection

A. N95 - Disposable filter (95% filter efficiency level) effective against particulate aerosols free of oil. Time use restrictions may apply.

2. Cautions and limitations

A. Not for use in atmospheres containing less than 19.5% oxygen.
B. Not for use in atmospheres immediately dangerous to life or health.
C. Do not exceed maximum use concentrations established by regulatory standards.
D. Failure to properly use and maintain this product could result in injury or death.
E. All approved respirators used for chemical, physical, and biological hazards must be maintained in accordance with MSDS, OSHA, and other applicable regulations.
F. Never wash, modify, melt, or cut parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
G. Refer to users instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
H. NIOSH will not evaluate respirators for use in surgical masks.
I. Special or critical user's instructions and/or specific use limitations apply. Refer to user's instructions (please download).

3. Special or critical user's instructions: The model RMK010 fitting box piece respirator has been manufactured by approved manufacturer for Company Kosto Products under TC 844-3004. Before use, the respirator must be fit tested either quantitatively or qualitatively by a qualified fit tester before using the respirator in a contaminated environment, as per the OSHA guidelines. (Fits box)

Produits Kosto / Kosto Products
60, Grande-Duvalville, Brossard, Québec, Canada J3C 6A3
1866-720-9292

NIOSH
N95
K
DECLASSIFIED

This respirator must be approved according to the following configurations:

TC	Filter/Cartridge	Model	Color and Markings
844-3004	N95	RMK010	X

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IMPORTANT NOTICE TO PURCHASER :
The following warranty is made in lieu of all other warranties expressed or implied, including but not limited to, the implied warranty of merchantability and fitness for purpose. Manufacturer warrants this product to be free from defects. Manufacturer's only obligation shall be, at its option, to replace any product proving defective or to refund the purchase price thereof. User shall determine the applicability of the product for its intended use, and assume all risks and liability in connection.

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USE AGAINST:
NIOSH approved for protection against particulate aerosols free of oil. Kosto recommends this respirator for use against particulate such as flour, iron dust and minerals.

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IMPORTANT LIMITATIONS: RESTRICTIONS IMPORTANTES :

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