USER'S INSTRUCTIONS

MYCO Medical Supplies, Inc.
N95 Respirator

MYTN01-11 (Medium)
MYTN01-12 (Small)

MYCO Medical Supplies, Inc.
2015 Production Drive, Apex, NC 27539
USA
TEL : +1-919-460-2535

TN01PL371Ulc
[Rev. c] 06/10/2020

USER'S INSTRUCTIONS

WARNING!

This respirator is intended to be worn to protect healthcare personnel from the transfer of microorganisms, body fluids and particulate material, however the risk concerning sickness or infection is not completely avoidable. Failure to read and follow all warnings and instruction will reduce respirator performance and may result in sickness and / or death. Before using, please read use instruction.

[USE INSTRUCTION]

This respirator does not supply oxygen and must not be used in atmospheres containing less than 19.5% oxygen.

Use for protection against particulate aerosols free of oil.

Do not use for gasses, vapors, asbestos, paint spray or any oil-based aerosols environments.

Before occupational use of this respirator a written respiratory protection program must be implemented meeting all the local government requirements. In the United States employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training, and fit testing.

Leave the use environment immediately in the case of dizziness, nausea, breathing difficulty or other symptoms of physical distress.

Don a new respirator if a proper fit cannot be obtained, if breathing becomes difficult, or the respirator becomes damaged. Only for single use.

Store unused respirators in a clean, cool, dark, dry place.

S-[SPECIAL OR CRITICAL USER'S INSTRUCTIONS]

These MYTN01-11 and MYTN01-12 filtering face piece respirators" have been manufactured by Approval Holder San-M Package Co., Ltd. for MYCO Medical Supplies, Inc. under TC-84A-3348.

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