

**1**

# FITTING INSTRUCTIONS

## ONE SIZE FITS ALL

### STEP 1

Hold the respirator in hand with the nosepiece at your fingertips, allowing the headbands to hang freely below your hand.

**2**

### STEP 2

Press the respirator firmly against your face with the nosepiece on the bridge of your nose.

**3**

### STEP 3

Stretch and position the top band high on the back of the head. Stretch the bottom band over the head and position below your ears.

**4**

### STEP 4

Using both hands, mold nosepiece to the shape of your nose.

**5**

### STEP 5

To test fit: **A)** Cup both hands over the respirator being careful not to disturb position, and **B)** inhale vigorously. If air leaks around the edges, reposition the straps or adjust strap tension for better fit. Please follow these fitting instructions during each use to achieve proper fit.