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ATTENTION

WARNING TO USER:

1. Follow all instructions and warnings on the use of this respirator and wear during all times of exposure. Failure to do so will reduce respirator effectiveness, wearer protection and may result in sickness or death.
2. Before use, the user must first be trained by the employer in proper respirator use, in accordance with applicable safety and health standards, for the contaminant and exposure level in the assigned work area.
3. The particulates and other contaminants which can be dangerous to your health include those which you cannot see, taste, or smell.
4. This product has not been sold with warnings or use instructions for personnel involved in health-care or related situations, where there may be the possibility of contact with disease or biological hazards. If you are considering such uses, first call the Moldex Technical Service Dept., +1 (310) 837-6500 or +1 (800) 421-0668, ext. 512/550.

RESTRICTIONS:

1. This respirator does not supply oxygen and must not be used in atmospheres containing less than 19.5% oxygen.
2. Do not remain in contaminated area if any physical distress occurs, for example breathing difficulty, dizziness or nausea.
3. Leave contaminated area and replace respirator if it is damaged, distorted, a proper fit can not be obtained or breathing becomes difficult.
4. Prior to each use, carefully inspect the entire respirator, including filter media and strap/weld attachment area for tears and damage. Staple or strap/weld attachment perforations do not affect respirator performance.
5. Do not alter, modify, or abuse this respirator.
6. Store unused respirators in box/bag in a clean, dry, non-contaminated area.
7. Dispose of respirator according to your employer's policy and local regulations.
8. Dispose no later than thirty days after first use.
9. Use respirator before the "use by" expiration date printed on box/bag.
10. If used for welding, wear appropriate eye and face protection.

For technical assistance call Moldex Technical Service Department.
+1 (310) 837-6500 or +1 (800) 421-0668, ext. 512/550.

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1. Particulate aerosols free of oil.

DO NOT USE AGAINST:

1. Concentrations of contaminants which are unknown, or are immediately dangerous to life or health.
2. Concentrations of particulates which exceed the maximum use concentration or 10 times the OSHA Permissible Exposure Limit (or Exposure Standard), whichever is lower.
3. Gases, vapors, asbestos, paint spray, sandblasting or particulate materials which generate harmful vapors.
4. Oil-based mists.

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Temperature range of storage conditions

<90%

Certified Product

Australian Standard

AS 1716 Lc 10876 S&H Global

In Australia Only, 4200 Approvals

This respirator has been manufactured to comply with the requirements of Australian Standards, AS 1716. Moldex-Metric, Inc. operates under Australian Standards guidelines for supervision, quality control and testing in accordance with the system of certification marking mandated by Australian Standards.

Approved as a Class P2 Particulate Respirator against mechanically and thermally generated dusts of low toxicity. Including free silica, coal, lead, iron ore, aluminum, manganese, talc, cotton, flour, and others found in mining, agriculture and manufacturing. For use against mists of materials with low toxicity which do not produce harmful gases or vapors. Also for use against metal fumes.

Contact your distributor or supplier for further information or phone Moldex U.S. at +1 (310) 837-6500 or fax +1 (310) 837-9563.

BIOHAZARD WARNINGS FOR MOLDEX PARTICULATE RESPIRATORS:

This respirator has not been approved for use against biohazards. You must consult with Australian Standards for regulations before use. Failure to follow these instructions could result in sickness or death.

LIMITATION: Respirators may reduce but not eliminate wearer exposure to airborne biohazards or risk of contracting any disease or infection. Only use this respirator as part of a comprehensive respiratory protection program. You will receive no respiratory protection if this respirator is not properly fitted and worn.

Moldex respirators must not be used on children.

For further information on use of respirator against biohazards contact Moldex USA at Fax +1 (310) 837-9563 ext. 512 or email to sales@moldex.com

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10 Respirators
10 Respiradores

FITTING

1. Users must follow instructions each time respirator is worn.
2. It's the wearer's responsibility to obtain a proper fit.
3. To choose the correct model/size Moldex Respirator you must fit test and user seal check.
4. If you cannot obtain a proper fit do not enter the contaminated area.
5. OSHA regulation 1910.134 (f) requires that the user must be fit tested before being assigned a respirator and annually thereafter.
6. Use Moldex Qualitative Fit Test Kit. User must be clean-shaven. Any facial hair, such as beards or long sideburns, may prevent the respirator from fitting properly.
7. Eyewear must not interfere with face to facepiece seal.
8. Re-fit test when there are any changes to the wearers physical condition that may affect fit.

1

2

3

4

5

1. Pull adjustment clip and strap fully below the bottom of the mask. (Fig 1)
2. Place strap around neck so outside of shell is against chest. (Fig 2)
3. a) Fit mask to face and pull top of strap to crown of head. Nose cushion must not be crossed inside of respirator.
b) Always use both top and bottom straps. (Fig 3)
4. To tighten pull both strap tabs. (Fig 4). To loosen, grasp both sides of neck strap and pull.
5. Each time user enters contaminated area respirator must be seal checked. Cover front of respirator by cupping both hands. (Fig 5) INHALE SHARPLY. A negative pressure should be felt inside respirator. If any leakage is detected at respirator edges, tighten or adjust strap by pulling back along the sides and/or reposition respirator. Repeat until sealed properly, otherwise see your supervisor. Entry into a contaminated area with an improper fit may result in sickness or death.
6. To hang around neck, pull top strap away from head, then pull mask down. (Fig 2)

CAUTION: When changing from any model/size to another respirator, you must fit test. If used against biohazards discard immediately after removing from face and do not hang around neck.

4601-802 REV A 01/12
6 x 5 x 8
CMYK

4601-802 REV A 01/12
6 x 5 x 8
CMYK