**Type MCRN951 NIOSH Approved Respirator**  
**Respirator Fitting Instructions:**  
**Instrucciones de Colocación del Respirador:**

1. Place fingers inside the respirator. Bend the nose clip around fingers to form nose shape.
2. Hold mask in position over nose and mouth with retaining straps around back of head. Pull lower portion of strap over head and fit around neck.
3. Stretch and pull the upper strap over your head on the back of your head above your ears.
4. Using both hands, form nose clip around nose and adjust the tension in the upper and lower portions of the continuous loop head strap, to ensure a positive face seal.
5. (a) To test fit for respirator without exhalation valve: Cup both hands over the respirator and exhale sharply. If air flows around your nose, tighten the nosepiece. If air leaks around the edges, reposition the headband for better fit.  
(b) To test fit for respirator with exhalation valve: Cup both hands over the respirators and inhale sharply. A negative pressure should be felt inside respirator. If air flows around your nose, tighten the nose-piece. If air leaks around the edges, reposition the headband for better fit.
6. Change respirator immediately if breathing becomes difficult or respirator becomes damaged or distorted, or a proper face fit cannot be maintained. Careful observance of these instructions is an important step in safe respirator use.

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