



FITTING INSTRUCTIONS:

Before fitting the respirator, ***carefully inspect the entire respirator*** to ensure that the respirator is not damaged and is functioning correctly. ***Discard the entire respirator if any part does not function properly.***

1. Place the respirator over your nose and mouth. Be sure the nose shape position is on top.
2. Pull the bottom strap over your head until it rests below the ears at the back of your neck. Untwist the strap.
3. Pull the top strap over your head until it rests on the crown of your head above your ears. Untwist the strap.
4. Without removing the respirator, adjust the straps until the respirator is secure. To tighten the respirator, gently pull the ends of both straps. To loosen the respirator, press on the back side of the buckle.
5. To check the correct fit, cup with both hands the respirator and exhale powerfully. If any airflow indicates leakage around the edges of the respirator stretch the headbands again by pulling them backwards. Then correct the position of the respirator (see picture) and check again. **WARNING!** It is not allowed to use the mask without proper fit. ***NEVER USE THE MASK AND DO NOT ENTER A CONTAMINATED AREA WITHOUT A PROPER SEAL.***

Note: Pinching the nosepiece using one hand may result in improper fit and less effective respirator performance. Use two hands.

Special Instructions: The model 3030V filtering facepiece respirator has been manufactured by Approval Holder Fido Masks CO., Ltd. for Alpha Solway Limited under TC-84A-8090.

Produce date: YYYY/MM

Lot No. XXXXXX

(the Lot No. XXXXXX will be as same as abbreviated Label indicated, the "XXXXXX" should be printed as purchase order number)

Contents: 10 pcs.

Fig. 4 Box printing side 3