

Fitting Instructions

Follow these instructions each time the mask is worn.

1. Pull open the top and bottom edges of the respirator to fully open. Grip head strap on either side of the slit in the middle of the strap and pull apart to create two separate head straps.
2. Bend the nose wire into a slight curve so that the respirator forms a cup.
3. Position the lower head strap around the back of the neck. Position the top head strap across the crown of the head.
4. Mold the nose area to the shape of the face, running the fingertips of both hands from the top of the nose-piece down both sides while pressing inward. Always use both hands.



Fit-check the Face-Seal as Follows:

1. Place both hands over the respirator - do not disturb the position of the respirator;
2. Exhale sharply (a positive pressure should be felt inside the respirator);
3. If you detect any air leaks around edges of the respirator, readjust the head straps and/or the nose-piece.

Warning:

If you cannot achieve satisfactory fit, do **NOT** enter the contaminated area. See your Supervisor.

User Instructions:

1. Failure to follow all instructions and limitations concerning the use of the respirator and/or failure to wear the respirator for the whole duration of exposure to contaminants can seriously reduce the performance characteristics of the respirator and lead to illness.
2. Before use, the wearer must first be trained by the employer in the correct use of the respirator in accordance with applicable safety and health standards.
3. OSHA Standards, including 29 CFR 1910.134(e)(5),(f) and CSA Standard Z94.4-93 require that the wearer be fit-tested.
4. Leave the contaminated area if breathing becomes difficult or if the respirator becomes damaged.
5. Discard the respirator and replace it with a new one if:
 - Excessive clogging of the respirator causes breathing difficulty
 - The respirator becomes damaged
 - The respirator becomes contaminated with blood or bodily fluids
6. Discard the respirator after every use.
7. Do **NOT** use the respirator with beards or other facial hair that interferes with direct contact between the face and the edge of the respirator, or any other conditions that may prevent a good face-seal.
8. This respirator is **NOT** to be used on children.

S- Special or Critical User's Instructions:

The model Actiprotect UF N95 Respirator has been manufactured by Sperian Respiratory Protection USA, LLC, Santa Ana, CA. for GlaxoSmithKline under TC-84A-5163.

Warning:

This respirator is intended to help reduce breathing in pathogenic biological airborne particulates or germs during seasonal influenza A or an influenza A pandemic. If you do not follow all instructions and limitations on the use of this respirator and/or do not wear this respirator during all times of exposure to the airborne germs the respirator may not be as effective. The respirator cannot stop the breathing in of all germs, in the air and does not eliminate the risk of disease or illness. Misuse may result in illness.

There is a high differential pressure associated with this mask. The air inside the respirator may feel hot when worn and may cause problems in individuals with underlying pulmonary disease.

Surgical masks, when worn properly, may reduce potential exposure of the wearer to blood and body fluids, but do not eliminate the risk of contracting any disease or infection.

FOR MORE INFORMATION or assistance call GlaxoSmithKline 1-800-825-5249.

Store below 30°C in original package