

Profit Royal 84A-7761

Use for:

Reduction of exposure to particles such as dust, pollens and particulate matter. Liquid or non-oil based particles from sprays that do not also emit oil aerosols or vapours

Do not use for:

Do not use for gases and vapours, oil aerosols, asbestos, or sandblasting; particulate concentrations that exceed either 10 times the occupational exposure limit or applicable government regulations, whichever is lower.

Warnings & Limitations

1) Before occupational use of this respirator a written respiratory protection program must be implemented meeting all the local government requirements. In the United States employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training, and fit testing.

2) The mask is designed to reduce exposures to certain airborne contaminants when fitted correctly. As such the mask is not suitable for kids, those with beards or other facial hair or other conditions that prevent a good seal between the face and the sealing surface of the respirator.

3) The Mask does not eliminate the risk of disease or infection nor should it be used as a substitute for proper health management. Wearers should still take the necessary precautions when entering and working in a hazardous setting. Failure to follow all instructions and limitations can reduce the masks effectiveness and may result in sickness or death.

4) The mask is designed as a single use disposable face mask - do not alter or wash this respirator.

5) Individuals with a compromised respiratory system, such as asthma or emphysema should consult a physician prior to use.

Time Use Limitation

If respirator becomes damaged, soiled, or breathing becomes difficult, leave the contaminated area immediately and replace the respirator.

Fitting



Storage Conditions :

Temperature | -20°C to +30°C
Humidity | <80% RH

存放環境

溫度 | -20°C 至 +30°C
濕度 | <80% 相對濕度

Note: The mask has an adjustable nose strip which is to be molded to conform to the user's nose. User should perform a User Seal Check prior to each wearing. To check the mask-to-face seal, place both hands completely over the mask and exhale. Be careful not to disturb the position of the mask. If air leaks around nose, readjust the nose clip as shown above. If air leaks at the mask edges, adjust the position of the mask. For more information contact the manufacturer below.



Manufacturer : Profit Royal Pharmaceutical Ltd
Website : www.nask.hk

Customer Services Hotline :
852-2387 8926 (Can) / 852-2387 8927 (Man) / 852-3153 4065 (Eng)